

A detailed description of the ANY-maze measures

ANY-maze
www.anymaze.com

Version 4, 30th May 2024

ANY-maze measures

1	Information measures.....	13
1.1	Test number	13
1.2	Animal number	13
1.3	Treatment	13
1.4	Treatment code.....	13
1.5	Stage.....	13
1.6	Trial number.....	14
1.7	Apparatus.....	14
1.8	The reason for test end.....	14
1.9	Test date	15
1.10	Test day of the week.....	15
1.11	Test time	15
1.12	Test time of day	15
1.13	User	16
1.14	Test notes.....	16
1.15	Animal notes	16
1.16	Animal lighter/darker than apparatus	16
1.17	Animal length	16
1.18	Percentage of frames tracked.....	17
1.19	Percentage of tracked frames head tracked.....	17
1.20	Source video file.....	17
1.21	Recorded video file	17
1.22	Video time when test started	18
1.23	Time period	18
1.24	Segment of test.....	18
1.25	Location of movable zones	18
2	Apparatus measures	19
2.1	Test duration.....	19
2.2	Time not hidden.....	19
2.3	Total distance travelled.....	19
2.4	Total distance travelled by the animal's head	20
2.5	First zone entered	20

2.6	Visited zone list	21
2.7	Investigated zone list	21
2.8	Average speed.....	21
2.9	Average speed when not hidden	22
2.10	Maximum speed	22
2.11	Total freezing episodes	22
2.12	Total time freezing.....	23
2.13	Latency to start of first freezing episode	23
2.14	Average freezing score.....	23
2.15	Total time mobile.....	23
2.16	Total time immobile.....	24
2.17	Total mobile episodes	24
2.18	Total immobile episodes.....	24
2.19	Latency to start of first mobile episode	25
2.20	Latency to start of first immobility episode.....	25
2.21	Latency to start of last mobile episode.....	25
2.22	Latency to start of last immobility episode.....	26
2.23	Total time active	26
2.24	Total time inactive	26
2.25	Total active episodes.....	27
2.26	Total inactive episodes.....	27
2.27	Longest active episode.....	27
2.28	Shortest active episode.....	28
2.29	Longest inactive episode.....	28
2.30	Shortest inactive episode.....	28
2.31	Rotations of the animal's body	29
2.32	Clockwise rotations of the animal's body	30
2.33	Anti-clockwise rotations of the animal's body.....	30
2.34	Partial rotations of the animal's body.....	30
2.35	Clockwise partial rotations of the animal's body.....	30
2.36	Anti-clockwise partial rotations of the animal's body	31
2.37	Absolute turn angle.....	31
2.38	Absolute head turn angle.....	31

2.39	Clockwise head turn angle	32
2.40	Anti-clockwise head turn angle.....	32
2.41	Number of rears.....	32
2.42	Total time rearing	32
2.43	Latency to first rear.....	33
2.44	Average duration of a rear.....	33
2.45	Maximum duration of a rear.....	33
2.46	Minimum duration of a rear	34
2.47	Path efficiency.....	34
2.48	Number of line crossings	34
2.49	On/off inputs positive reversal	35
2.50	On/off inputs negative reversal	35
2.51	Tracking quality.....	36
2.52	Number of centre positions recorded	36
2.53	Number of head positions recorded.....	37
2.54	RAPC - Type 1 errors	37
2.55	RAPC - Type 2 errors	37
2.56	RAPC - Door sequence	38
3	Zone measures.....	39
3.1	Number of entries to the zone	39
3.2	Number of exits from the zone.....	39
3.3	Number of entries of the animal's head into the zone.....	39
3.4	Time in the zone.....	39
3.5	Time the animal's head was in the zone.....	40
3.6	Was first zone entered.....	40
3.7	Number of bouts of investigation of the zone.....	40
3.8	Total time investigating the zone	41
3.9	Was first zone investigated.....	41
3.10	Longest bout of investigation of the zone	41
3.11	Shortest bout of investigation of the zone	42
3.12	Average duration of investigation of the zone	42
3.13	List of the duration of each investigation of the zone	42
3.14	Number of partial exits from the zone	43

3.15	Time partially exited from the zone.....	43
3.16	Distance travelled in the zone	44
3.17	Distance travelled until first entry into the zone	45
3.18	Distance travelled by the animal's head in the zone	45
3.19	Distance travelled while investigating the zone	46
3.20	Distance travelled before first investigation of the zone	46
3.21	Latency to first entry to the zone.....	46
3.22	Latency to first exit from the zone.....	47
3.23	Latency to last entry to the zone	47
3.24	Latency to the first entry of the animal's head into the zone	48
3.25	Latency to the first exit of the animal's head from the zone.....	48
3.26	Latency to first investigation of the zone	48
3.27	Latency to end of first investigation of the zone	49
3.28	Average speed in the zone.....	49
3.29	Maximum speed in the zone.....	49
3.30	Average speed while investigating the zone.....	50
3.31	Longest visit to the zone	50
3.32	Shortest visit to the zone	51
3.33	Average duration of visit to the zone	51
3.34	List of the duration of each visit to the zone	51
3.35	Time mobile in the zone	52
3.36	Time immobile in the zone	52
3.37	Immobile episodes in the zone	52
3.38	Time mobile while investigating the zone	53
3.39	Time immobile while investigating the zone	53
3.40	Immobile episodes while investigating the zone.....	54
3.41	Time active in the zone	54
3.42	Time inactive in the zone	55
3.43	Inactive episodes in the zone.....	55
3.44	Initial distance from the zone	56
3.45	Average distance from the zone	56
3.46	Maximum distance from the zone.....	57
3.47	Minimum distance from the zone	57

3.48	Cumulative distance from the zone	58
3.49	Average distance of the animal's head from the zone	58
3.50	Maximum distance of the animal's head from the zone	59
3.51	Minimum distance of the animal's head from the zone.....	59
3.52	Average distance to the zone border	59
3.53	Maximum distance to the zone border	60
3.54	Minimum distance to the zone border	61
3.55	Average distance from the animal's head to the zone border	61
3.56	Maximum distance from the animal's head to the zone border	62
3.57	Minimum distance from the animal's head to the zone border.....	62
3.58	Time getting closer to the zone	62
3.59	Time getting further away from the zone.....	63
3.60	Initial heading error to the zone	64
3.61	Signed initial heading error to the zone	66
3.62	Average absolute heading error to the zone	67
3.63	Time moving towards the zone	69
3.64	Time moving away from the zone	71
3.65	Time oriented towards the zone.....	72
3.66	Time oriented towards the centre of the zone when inside zone.....	73
3.67	Absolute turn angle while in the zone	73
3.68	Absolute head turn angle while in the zone	74
3.69	Freezing bouts in the zone.....	74
3.70	Time freezing in the zone.....	75
3.71	Freezing bouts while investigating the zone.....	75
3.72	Time freezing while investigating the zone	75
3.73	Number of rears in the zone	75
3.74	Total time rearing in the zone.....	76
3.75	Latency to first rear in the zone	76
3.76	Average duration of a rear in the zone	76
3.77	Maximum duration of a rear in the zone.....	77
3.78	Minimum duration of a rear in the zone	77
3.79	Time spent in Whishaw's Corridor	77
3.80	Distance travelled in Whishaw's Corridor.....	78

3.81	Time the animal's head was in the zone when its centre was outside the zone.....	79
3.82	Path efficiency to first entry to the zone	80
3.83	Corrected integrated path length	80
3.84	Number of line crossings while in the zone	81
4	Point measures	82
4.1	Average distance from the point	82
4.2	Maximum distance from the point	82
4.3	Minimum distance from the point.....	82
4.4	Average distance of the animal's head from the point	83
4.5	Maximum distance of the animal's head from the point	83
4.6	Minimum distance of the animal's head from the point	83
4.7	Time moving towards the point.....	83
4.8	Time moving away from the point.....	84
4.9	Time the animal's head was moving towards the point	85
4.10	Time the animal's head was moving away from the point	86
4.11	Average speed moving towards the point.....	87
4.12	Initial heading error to the point	88
4.13	Average absolute heading error to the point	89
4.14	Time the animal's head was oriented towards the point	89
4.15	Time the animal's head was oriented away from the point	90
4.16	Number of times the animal's head was oriented towards the point.....	91
4.17	Point X coordinate	92
4.18	Point Y coordinate.....	93
4.19	Approximate time at the point	93
5	Sequence measures	94
5.1	Number of sequences	94
5.2	Total time performing sequences	94
5.3	Latency to start of first sequence	94
5.4	Latency to completion of first sequence	95
5.5	Average time to complete the sequence.....	95
5.6	Maximum duration of a sequence.....	95
5.7	Minimum duration of a sequence	95
5.8	Total distance travelled during sequences	96

5.9	Average distance travelled per sequence.....	96
5.10	Maximum distance travelled during a sequence.....	96
5.11	Minimum distance travelled during a sequence	97
5.12	Average speed during the sequence.....	97
6	Key measures	98
6.1	Number of presses	98
6.2	Time pressed	98
6.3	Latency to first press	98
6.4	Latency to first release.....	99
6.5	Distance travelled before 1st press	99
6.6	Longest press	99
6.7	Shortest press	100
6.8	Average press duration.....	100
6.9	Press frequency.....	100
6.10	List of activation durations	101
7	On/off input measures.....	102
7.1	Number of activations	102
7.2	Time active.....	102
7.3	Latency to first activation	102
7.4	Latency to first deactivation	103
7.5	Longest activation.....	103
7.6	Shortest activation.....	103
7.7	Average activation duration.....	104
7.8	Frequency of activations.....	104
8	Rotary encoder measures.....	105
8.1	Number of rotations	105
8.2	Time turning.....	105
8.3	Distance (only applies to running wheels).....	105
8.4	Number of clockwise rotations.....	106
8.5	Number of anti-clockwise rotations	106
8.6	Number of reversals.....	106
8.7	Number of half rotations	106
8.8	Number of quarter rotations	107

8.9	Degrees of clockwise rotation.....	107
8.10	Degrees of anti-clockwise rotation	108
8.11	Maximum RPM.....	108
8.12	Minimum RPM	109
8.13	Average RPM.....	109
8.14	Average RPM while turning	109
9	Signal measures	111
9.1	Average	111
9.2	Maximum	111
9.3	Minimum.....	111
9.4	Time of maximum	111
9.5	Time of minimum.....	112
9.6	Baseline.....	112
9.7	Baseline deviation.....	112
9.8	Baseline standard deviation.....	113
9.9	Time of the end of the baseline period.....	113
9.10	Time of first positive deviation from baseline	113
9.11	Time of first return to baseline from positive deviation.....	113
9.12	Time of first negative deviation from baseline	114
9.13	Time of first return to baseline from negative deviation	114
9.14	Integral above baseline.....	114
9.15	Integral below baseline.....	115
9.16	Average maximum for each visit to zone.....	115
9.17	Average time to maximum for each visit to zone	115
9.18	Average minimum for each visit to zone	116
9.19	Average time to minimum for each visit to zone.....	116
9.20	Average at zone entry	117
9.21	Average at zone exit.....	117
10	Sensor measures	118
10.1	Initial value.....	118
10.2	Average value.....	118
10.3	Maximum value	118
10.4	Minimum value	118

10.5	Change	119
11	Movement detector measures	120
11.1	Count.....	120
11.2	Time moving.....	120
11.3	Latency to first beam break	120
12	On/off output measures	122
12.1	Number of activations	122
12.2	Time active	122
12.3	Latency to first activation	122
12.4	Latency to first deactivation	123
12.5	Longest activation	123
12.6	Shortest activation	123
12.7	Average activation duration.....	123
12.8	Frequency of activations.....	124
13	Speaker measures	125
13.1	Number of activations	125
13.2	Time active.....	125
13.3	Latency to first activation	125
13.4	Latency to first deactivation	126
13.5	Longest activation	126
13.6	Shortest activation	126
13.7	Average activation duration.....	127
13.8	Frequency of activations.....	127
14	Shocker measures	128
14.1	Number of activations	128
14.2	Time active.....	128
14.3	Latency to first activation	128
14.4	Latency to first deactivation	129
14.5	Longest activation	129
14.6	Shortest activation	129
14.7	Average activation duration.....	130
14.8	Frequency of activations.....	130
15	Pellet dispenser measures	131

15.1	Number of pellets dispensed	131
15.2	Latency to first pellet dispensed	131
16	Laser controller measures.....	132
16.1	Number of activations	132
16.2	Time active.....	132
16.3	Latency to first activation	132
16.4	Latency to first deactivation	133
16.5	Longest activation.....	133
16.6	Shortest activation.....	133
16.7	Average activation duration.....	134
16.8	Frequency of activations.....	134
17	Syringe pump measures.....	135
17.1	Volume infused	135
17.2	Volume withdrawn.....	135
18	Light controller measures	136
18.1	Number of activations	136
18.2	Time active.....	136
18.3	Latency to first activation	136
18.4	Latency to first deactivation	137
18.5	Longest activation.....	137
18.6	Shortest activation.....	137
18.7	Average activation duration.....	138
18.8	Frequency of activations.....	138
19	OPAD measures	139
19.1	Temperature when contact broken.....	139
19.2	Number of non-lick contacts.....	139
19.3	Temperature of interest: time in contact	139
19.4	Temperature of interest: number of times contact broken	140
19.5	Temperature of interest: number of times contact made	140
19.6	Temperature of interest: number of licks.....	140
20	Procedure measures	141
20.1	Value	141
20.2	Average value.....	141

20.3	Maximum value	141
20.4	Minimum value	142
20.5	Sum of values	142
20.6	Count of values	143
20.7	List of values.....	143
21	Event measures.....	144
21.1	Number of events	144
21.2	Latency to first event	144
22	Virtual switch measures.....	145
22.1	Number of activations	145
22.2	Time active	145
22.3	Latency to first activation	145
22.4	Latency to first deactivation	146
22.5	Distance travelled before 1st activation	146
22.6	Distance travelled while active	146
22.7	Longest activation	147
22.8	Shortest activation	147
22.9	Average activation duration.....	147
22.10	Frequency of activations.....	148

1 Information measures

ANY-maze includes the following information measures for each test:

1.1 Test number

<i>Description</i>	A sequential number assigned to each test when it's completed.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.2 Animal number

<i>Description</i>	The animal number of the animal this test was performed on.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	Animal numbers are assigned sequentially to animals as they're added to an experiment. These numbers are never reissued, even if an animal is removed from the experiment. For example, if you create an experiment with 10 animals in it, they'll be numbered 1-10. If you immediately delete the 10 animals and then add 10 new animals, the new ones will be numbered 11-20 and the experiment won't have any animals numbered 1-10 at all.

1.3 Treatment

<i>Description</i>	The name of the treatment given to the animal this test was performed on.
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.4 Treatment code

<i>Description</i>	The code of the treatment given to the animal this test was performed on.
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.5 Stage

<i>Description</i>	The name of the stage which this test was part of.
<i>Independent variable</i>	Yes

<i>Dependent variable</i>	No
<i>Notes</i>	All tests must be part of a stage - even if an experiment only contains a single default stage.





1.6 Trial number

<i>Description</i>	The trial number within the stage for the animal which this test was performed on.
<i>Independent variable</i>	Yes - see notes.
<i>Dependent variable</i>	No
<i>Notes</i>	<p>Within a stage, an animal can be tested more than once - i.e. have repeated trials. The first test in a stage is the animal's trial 1, the second test is the animal's trial 2, etc.</p> <p>Although the trial numbers will be repeated in different stages, ANY-maze still views them as different when using Trial number as an independent variable. Thus, for example, Trial 1 in an <i>Acquisition</i> stage would be seen as different to trial 1 in a <i>Re-test</i> stage.</p>

1.7 Apparatus

<i>Description</i>	The name of the apparatus that the test was performed on.
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.8 The reason for test end

<i>Description</i>	<p>The reason the test ended. This is either a description of a standard reason (see notes), the name of a procedure that ended the test, or a test end reason created for a procedure.</p> <p> <i>If your protocol is set up to use Events and Actions rather than procedures, then the test end reason can be the name of the event that triggered an action that ended the test (the name of the event is used because this reflects the reason why the test ended, for example, 'Island entry').</i></p>						
<i>Independent variable</i>	Yes						
<i>Dependent variable</i>	Yes. Nominal data. Can't be analysed across time.						
<i>Notes</i>	<p>The standard reasons for ANY-maze to end a test are:</p> <table> <tr> <td><i>User ended test</i></td> <td>The user ended the test by clicking the  <i>Stop test</i> button.</td> </tr> <tr> <td><i>Test duration</i></td> <td>The test ended because the test duration was reached.</td> </tr> <tr> <td><i>End of video</i></td> <td>This reason will be used if you're using a video file for the test (rather than tracking the animal 'live'), and the end of this video file is reached.</td> </tr> </table>	<i>User ended test</i>	The user ended the test by clicking the  <i>Stop test</i> button.	<i>Test duration</i>	The test ended because the test duration was reached.	<i>End of video</i>	This reason will be used if you're using a video file for the test (rather than tracking the animal 'live'), and the end of this video file is reached.
<i>User ended test</i>	The user ended the test by clicking the  <i>Stop test</i> button.						
<i>Test duration</i>	The test ended because the test duration was reached.						
<i>End of video</i>	This reason will be used if you're using a video file for the test (rather than tracking the animal 'live'), and the end of this video file is reached.						

<i>Lost video signal</i>	The video signal showing the apparatus was lost (for example, the camera was unplugged).
<i>Lost I/O device</i>	An I/O device used in the experiment was lost (i.e. unplugged from the ANY-maze computer).
<i>Out of memory</i>	While running a test, ANY-maze records the test results in the computer's memory - if this memory becomes full, then the test will be ended for this reason. This is very unlikely to occur.
<i>Error saving results</i>	For some reason, ANY-maze is unable to save the results to the experiment file. This is very unlikely to occur.

1.9 Test date

<i>Description</i>	The date when the test was <i>started</i> .
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.10 Test day of the week

<i>Description</i>	The day of the week (Monday-Sunday) when the test <i>started</i> .
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	Yes. Nominal data. Can't be analysed across time.
<i>Notes</i>	None

1.11 Test time

<i>Description</i>	The time when the test <i>started</i> .
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.12 Test time of day

<i>Description</i>	The time of day, am or pm, when the test <i>started</i> .
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	Yes. Nominal data. Can't be analysed across time.
<i>Notes</i>	None

1.13 User

<i>Description</i>	The name of the user who was logged on while the test was performed.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	None

1.14 Test notes

<i>Description</i>	Any notes recorded for the individual test - see notes.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	Only the first 80 characters of the notes are actually used for this measure's value. The notes themselves can be up to 32,000 characters in length.

1.15 Animal notes

<i>Description</i>	Any notes recorded for the individual animal this test was performed on - see notes.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	Only the first 80 characters of the notes are actually used for this measure's value. The notes themselves can be up to 32,000 characters in length.

1.16 Animal lighter/darker than apparatus

<i>Description</i>	Whether the animal is lighter or darker than the apparatus background, as specified in the protocol list under Animal colour.
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	Yes
<i>Notes</i>	None

1.17 Animal length

<i>Description</i>	The length of the animal in millimetres
<i>Independent variable</i>	No
<i>Dependent variable</i>	Yes
<i>Notes</i>	The animal length measure is only available if the animal length is required to be <i>entered*</i> and the protocol specifies that the animals have different lengths. *This is only the case when the protocol includes at least one hidden zone, or freezing is being detected but the animal is not being tracked.

1.18 Percentage of frames tracked

<i>Description</i>	The percentage of the frames (video pictures) sent by the camera (or video file) in which ANY-maze successfully tracked the animal.
<i>Independent variable</i>	No
<i>Dependent variable</i>	Yes
<i>Notes</i>	ANY-maze may not be able to determine the position of the animal in all the pictures it receives from the camera (or video file). For example, there might be a moving reflection (common in the water-maze) which interferes with tracking, or the animal might move into a dead space (where the camera can't fully see it), etc.

1.19 Percentage of tracked frames head tracked

<i>Description</i>	The percentage of the successfully tracked frames (video pictures) in which ANY-maze tracked the animal's head.
<i>Independent variable</i>	No
<i>Dependent variable</i>	Yes
<i>Notes</i>	ANY-maze will only try to detect the animal's head in frames in which the animal itself has been tracked, but it may not always be successful. The head tracking works by analysing the shape of the animal and if it is in a position which makes it look quite un-rodent like (for example, it has hunched up into a ball, or it is partially obscured by an object in a NOR test, etc.) then the head tracking may not be able to determine where the animal's head is.

1.20 Source video file

<i>Description</i>	The name of the video file from which a test was run. The full path will be displayed on the Data page; the file name only will be shown on the test schedule report.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	If the test was not run from a video, this will not be displayed.

1.21 Recorded video file

<i>Description</i>	The name of the video file that was recorded while a test was being run. The full path will be displayed on the Data page; the file name only will be shown on the test schedule report.
<i>Independent variable</i>	No
<i>Dependent variable</i>	No
<i>Notes</i>	If no video was recorded of this test, this will not be displayed.


1.22 Video time when test started

<i>Description</i>	The time in the video file at which the test started.
<i>Independent variable</i>	No
<i>Dependent variable</i>	Yes
<i>Units</i>	Seconds
<i>Notes</i>	This value will be reported as undefined for tests that were not performed from a video, for example tests performed from a live picture from a camera. It will also be reported as undefined for tests performed in versions prior to 6.33 of ANY-maze, as this was the first version to record this information.

1.23 Time period

<i>Description</i>	A Time period of the test
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	For more information about Time periods refer to An introduction to Time periods

1.24 Segment of test

<i>Description</i>	A segment of the test, described as x-y seconds; for example, 30-60 seconds.
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	No
<i>Notes</i>	<p>ANY-maze can break tests down into equal length segments and then perform analysis on them. The length of the segments that tests are broken into is specified in the protocol's Analysis across time element.</p> <p> It's very important to understand that when you use the <i>segment of test</i> measure, you effectively tell ANY-maze to start analysing segments rather than tests. This means, for example, that the Data page spreadsheet will include one row per test segment rather than one row per test.</p>

1.25 Location of movable zones

<i>Description</i>	The location of a movable zone in the test
<i>Independent variable</i>	Yes
<i>Dependent variable</i>	Yes
<i>Notes</i>	There is one measure for each movable zone in the protocol. You'll find more information about movable zone here.

2 Apparatus measures

2.1 Test duration

<i>Description</i>	Reports the duration of a test.
<i>Calculation method</i>	The value of the test clock when the test ended.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the amount of the test duration which fell inside the period. This will be the full duration of the period for all periods except the last one in the test. This result is most useful in calculations.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.2 Time not hidden

<i>Description</i>	Reports the amount of time in the test for which the animal was not hidden.
<i>Calculation method</i>	Each time the animal becomes hidden the amount of time it has been visible for is added to a running total of the time not hidden.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the amount of time during the period for which the animal was not hidden.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the protocol includes at least one hidden zone.

2.3 Total distance travelled

<i>Description</i>	Reports the total distance that the animal travelled during the test.
<i>Calculation method</i>	Sum of the distance between each point in the track - see note below.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the distance that the animal travelled during the period.
<i>Units</i>	Metres
<i>Notes</i>	In some situations, tracks can have small oscillations in them which tend to generate unrepresentatively large values for distance travelled. This occurs most often when an animal travels slowly while moving its body a lot - for example, while exploring an open field. To overcome this, ANY-maze uses an adaptive smoothing algorithm to attenuate these oscillations when calculating distance travelled - see figure 1. Note: The definition of what's a <i>small</i> oscillation is based on the animal's size.




Figure 1. Measuring the length of the actual track (shown in orange) would yield an unrepresentatively large value for distance travelled. ANY-maze uses a 'smoothed' track (shown in green) to better estimate the true distance travelled. [Note: The oscillations in this track have been exaggerated to aid explanation.]

This measure only includes the distance the animal travelled while ANY-maze could see it. In other words, if the animal was hidden for part of the test then the distance it may have travelled while hidden is not included in this value.

2.4 Total distance travelled by the animal's head

<i>Description</i>	Reports the total distance that the animal's head travelled during the test.
<i>Calculation method</i>	Sum of the distance between each point in the head track.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	As for the Total distance travelled, ANY-maze will smooth the animal's head track to remove small oscillations which would otherwise distort the result of this measure.

2.5 First zone entered

<i>Description</i>	Reports the name of the first zone that the animal entered during the test.
<i>Calculation method</i>	Simply reflects the first zone entry. This is affected by the <i>Don't score any results in this zone until the first 'true' entry</i> option on the <i>Zone entry settings</i> page. See <i>Choosing how ANY-maze should detect entries into a zone</i> for more details.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	When analysed, this value will be treated as a nominal value - see statistical tests included in ANY-maze. In some circumstances, it's possible that two or more zones could be the 'first zone entered'. For example, if you create a protocol in which a single area of the

apparatus is included in two zones then, if the animal enters this area first, it will have entered both zones at the same time - meaning there are two 'first zones entered'. In this situation, ANY-maze will report the 'first zone entered' as being the first one of the zones in the zone list shown in the protocol.

This measure is only available if no zones are excluded from analysis.

2.6 Visited zone list

<i>Description</i>	Reports a comma-separated list of the names of the zones the animal visited, in the order in which they were visited.
<i>Calculation method</i>	Each time the animal enters a zone, the zone's name is added to the list.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	It is possible to define an area of the apparatus as being part of two (or more) zones. If the animal enters such an area, then it is necessarily entering all the zones simultaneously and the list will therefore include them all. In this case, the zones are added to the list in the order in which they appear in the protocol. This measure is only available if no zones are excluded from analysis.

2.7 Investigated zone list

<i>Description</i>	Reports a comma-separated list of the names of the zones the animal investigated, in the order in which they were investigated.
<i>Calculation method</i>	Each time the animal investigates a zone, the zone's name is added to the list.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	It is possible for an animal to start investigating two (or more) zones at exactly the same time. If this occurs then the zones are added to the list in the order in which they appear in the protocol. This measure is only available if no zones are excluded from analysis.

2.8 Average speed

<i>Description</i>	Reports the average speed of the animal during a test.
<i>Calculation method</i>	Calculated by dividing the Total distance travelled by the Test duration.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the Total distance travelled during the period divided by the Test duration for the period. Note: The 'Test duration' for a period is the amount of the test duration which fell in the period - this is the period's duration for all periods except the last one in the test.
<i>Units</i>	Metres per second

Notes If you want to know average speed while mobile (i.e. ignoring periods when the animal was stationary), then use a calculation of *Total distance travelled / Total time mobile*.

2.9 Average speed when not hidden

Description Reports the average speed of the animal for the periods of the test for which it was not hidden.

Calculation method Calculated by dividing the Total distance travelled by the Time not hidden.

Analysis across time This measure can be analysed across time. The result for a period is the Total distance travelled during the period divided by the Time not hidden for the period. Note: If the Time not hidden for a period is zero (i.e. the animal was hidden for the entire period) then the result for this measure will be undefined.

Units Metres per second

Notes This measure is only available if the protocol includes at least one hidden zone.

2.10 Maximum speed

Description Reports the maximum speed of the animal.

Calculation method The speed of the animal between positions is calculated and the maximum speed is found.

Analysis across time This measure can be analysed across time.

Units Metres per second

Notes The calculation of maximum speed does not use *successive* positions but instead requires that the animal move at least a minimum distance (which is based on the animal's size) and the speed to cover *this* distance is calculated. This method of calculation is used to avoid reporting the speed of movements that don't constitute locomotion of the animal. For example, if an animal scratches, its centre point may oscillate rapidly but this will not be reported as the animal's maximum speed.

2.11 Total freezing episodes

Description Reports the number of times the animal froze during the test.

Calculation method Each time the animal begins to freeze, a counter is incremented. The result is the counter's value at the end of the test.

Analysis across time This measure can be analysed across time. The result for a time period is the number of freezing episodes that *started* during the period.

Units None

Notes None

2.12 Total time freezing

<i>Description</i>	Reports the total amount of time during the test that the animal was freezing.
<i>Calculation method</i>	The duration of each freezing episode is calculated and these values are summed.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the amount of time the animal was freezing during the period. If the animal is freezing at the start of the period, then the result includes the time until the animal stops freezing (or the period ends). This means that the total freezing episodes during a time period can be zero when the time freezing during the period is non-zero.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.13 Latency to start of first freezing episode

<i>Description</i>	Reports the latency to the start of the first moment in the test when the animal freezes. If the animal is freezing at the start of the test, this value will be zero.
<i>Calculation method</i>	The test time when the first freezing episode occurs.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the time when the animal <i>first</i> froze during the time period. This means that if the animal is freezing at the start of the time period, then the latency is not reported as zero.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.14 Average freezing score

<i>Description</i>	Reports the average of the animal's freezing score.
<i>Calculation method</i>	The freezing scores are summed and divided by their count.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	This measure is most useful when applied to time segments, as comparison of the score between segments gives an indication of the changes in the animal's activity across time.

2.15 Total time mobile

<i>Description</i>	Reports the amount of time the animal was mobile during the test.
<i>Calculation method</i>	Calculated by subtracting the Total time immobile from the Test duration.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Test duration of the period minus Total time immobile during the period. Note: The 'Test duration' for a period is the amount of the test duration which fell in the period - this is the period's duration for all periods except the last one in the test.
<i>Units</i>	Seconds

Notes None

2.16 Total time immobile

Description Reports the amount of time the animal was immobile during the test.

Calculation method Sums the duration of each immobile episode in the test. The definition of immobility depends on the current tracking options - see Detecting immobility.

Analysis across time This measure can be analysed across time. For any time period, the result is the sum of the duration of each immobile episode in the period.

Episodes of immobility which fall partly in a period but which start or end outside it are calculated as if they started or ended at the start or end of the period, respectively. This means that it is possible to have a result for Total immobile episodes in the period which is zero and a result for *Total time immobile* in the period which is not zero. For example, if the animal is immobile at the start of a period and remains immobile throughout the period then the *Total time immobile* in the period will be the period's duration but the Total immobile episodes in the period will be zero because no transition from a mobile state to an immobile state occurred during the period.

Units Seconds

Notes None

2.17 Total mobile episodes

Description Reports the number of times the animal was mobile during the test.

Calculation method Counts the number of transitions from an immobile state to a mobile state during the test. For the purposes of this calculation, the animal is assumed to be immobile at the start of the test. The definition of immobility depends on the current tracking options - see Detecting immobility.

Analysis across time This measure can be analysed across time. For any time period, the result is the count of the number of transitions from an immobile state to a mobile state during the period.

Units Seconds

Notes None

2.18 Total immobile episodes

Description Reports the number of times the animal was immobile during the test.

Calculation method Counts the number of transitions from a mobile state to an immobile state during the test. For the purposes of this calculation, the animal is assumed to be mobile at the start of the test. The definition of immobility depends on the current tracking options - see Detecting immobility.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the count of the number of transitions from a mobile state to an immobile state during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.19 Latency to start of first mobile episode

<i>Description</i>	Reports the latency to the start of the first moment in the test when the animal is considered to be mobile. If the animal is mobile at the start of the test, this measure's value will be zero.
<i>Calculation method</i>	The test time when the animal first becomes mobile.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.20 Latency to start of first immobility episode

<i>Description</i>	Reports the latency to the start of the first moment in the test when the animal is considered to be immobile. If the animal is immobile at the start of the test, this measure's value will be zero.
<i>Calculation method</i>	The test time when the animal first becomes immobile.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.21 Latency to start of last mobile episode

<i>Description</i>	Reports the latency to the start of the last episode of mobility in the test. If the animal is only mobile once in the test, this value will be the same as the <i>Latency to start of first mobile episode</i> . If the animal is mobile throughout the test then this value will be zero. If the animal is never mobile during the test then this value will be undefined.
<i>Calculation method</i>	The test time when the animal becomes mobile for the last time in the test.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.22 Latency to start of last immobility episode

<i>Description</i>	Reports the latency to the start of the last episode of immobility in the test. If the animal is only immobile once in the test, this value will be the same as the <i>Latency to start of first immobile episode</i> . If the animal is immobile throughout the test then this value will be zero. If the animal is never immobile during the test then this value will be undefined.
<i>Calculation method</i>	The test time when the animal becomes immobile for the last time in the test.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

2.23 Total time active

<i>Description</i>	Reports the amount of time the animal was active during the test.
<i>Calculation method</i>	Calculated by subtracting the Total time inactive from the Test duration.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Test duration of the period minus Total time inactive during the period. Note: The 'Test duration' for a period is the amount of the test duration which fell in the period - this is the period's duration for all periods except the last one in the test.
<i>Units</i>	Seconds
<i>Notes</i>	An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.24 Total time inactive

<i>Description</i>	Reports the amount of time the animal was inactive during the test.
<i>Calculation method</i>	Sums the duration of each inactive episode in the test.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the sum of the duration of the inactive episodes in the period. Episodes of inactivity which fall partly in a period, but which start or end outside it, are calculated as if they started or ended at the start or end of the period, respectively. This means that it is possible to have a result for Total inactive episodes in the period which is zero and a result for <i>Total time inactive</i> in the period which is not zero. For example, if the animal is inactive at the start of a period and remains inactive throughout the period, then the <i>Total time inactive</i> in the period will be the period's duration, but the Total inactive episodes in the period will be zero (because no transition from an active state to an inactive state occurred during the period).
<i>Units</i>	Seconds
<i>Notes</i>	Inactivity is defined as NOT activity. An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an

activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.25 Total active episodes

<i>Description</i>	Reports the number of times the animal was active during the test.
<i>Calculation method</i>	Counts the number of transitions from an inactive to an active state. For the purposes of this calculation, the animal is assumed to be inactive at the start of the test.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of transitions from an inactive to an active state that occurred during the period.
<i>Units</i>	None
<i>Notes</i>	An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.26 Total inactive episodes

<i>Description</i>	Reports the number of times the animal was inactive during the test.
<i>Calculation method</i>	Counts the number of transitions from an active to an inactive state. For the purposes of this calculation, the animal is assumed to be active at the start of the test.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of transitions from an active to an inactive state that occurred during the period.
<i>Units</i>	None
<i>Notes</i>	Inactivity is defined as NOT activity. An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.27 Longest active episode

<i>Description</i>	Reports the duration of the longest continuous period of activity during the test.
<i>Calculation method</i>	The duration of each episode of activity is calculated when the episode ends. The longest one is found.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the longest episode of activity that occurred during the period.

<i>Units</i>	Seconds
<i>Notes</i>	An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.28 Shortest active episode

<i>Description</i>	Reports the duration of the shortest continuous period of activity during the test.
<i>Calculation method</i>	The duration of each episode of activity is calculated when the episode ends. The shortest one is found.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest episode of activity that occurred during the period.
<i>Units</i>	Seconds
<i>Notes</i>	An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.29 Longest inactive episode

<i>Description</i>	Reports the duration of the longest continuous period of inactivity during the test.
<i>Calculation method</i>	The duration of each episode of inactivity is calculated when the episode ends. The longest one is found.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the longest episode of inactivity that occurred during the period.
<i>Units</i>	Seconds
<i>Notes</i>	Inactivity is defined as NOT activity. An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.30 Shortest inactive episode

<i>Description</i>	Reports the duration of the shortest continuous period of inactivity during the test.
<i>Calculation method</i>	The duration of each episode of inactivity is calculated when the episode ends. The shortest one is found.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest episode of inactivity that occurred during the period.
<i>Units</i>	Seconds

Notes

Inactivity is defined as NOT activity. An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the immobility detection element specifies that mobility should NOT be detected, then activity analysis will be based purely on the performance of other behaviours.

2.31 Rotations of the animal's body

Description Reports the number of times the animal's body completed an entire rotation of 360°.

Calculation method The animal's centre point is taken as a virtual origin; i.e. this origin is adjusted to be in the same place in each frame. A line is then taken from the animal's centre point to its head creating a vector. The angle between successive vectors is calculated and while the angle continues to have the same sign, the angles are accumulated - when the accumulated angle reaches 360°, the animal has completed a rotation. In fact this is a simplification, as the exact method used takes partial reversals of direction into account - see the figure below.




Figure 2. A partial reversal in direction doesn't alter the end of the rotation. A complete reversal however, means the animal has to rotate back to the reversal point to complete the rotation.

Analysis across time This measure can be analysed across time. A rotation is deemed to occur at the time it is completed.

Units None

Notes None

2.32 Clockwise rotations of the animal's body

<i>Description</i>	Reports the number of times the animal's body completed an entire rotation of 360° in a clockwise direction.
<i>Calculation method</i>	This measure is calculated in the same way as rotations of the animal's body only it reports just the clockwise rotations.
<i>Analysis across time</i>	This measure can be analysed across time. A rotation is deemed to occur at the time it is completed.
<i>Units</i>	None
<i>Notes</i>	None

2.33 Anti-clockwise rotations of the animal's body

<i>Description</i>	Reports the number of times the animal's body completed an entire rotation of 360° in an anti-clockwise direction.
<i>Calculation method</i>	This measure is calculated in the same way as rotations of the animal's body, except that it reports only the anti-clockwise rotations.
<i>Analysis across time</i>	This measure can be analysed across time. A rotation is deemed to occur at the time it is completed.
<i>Units</i>	None
<i>Notes</i>	None

2.34 Partial rotations of the animal's body

<i>Description</i>	Reports the number of times the animal's body rotated by at least the 'Partial rotation angle' defined in the protocol but didn't complete an entire rotation.
<i>Calculation method</i>	*** Feature still in beta ***
<i>Analysis across time</i>	This measure can be analysed across time. A partial rotation is deemed to occur at the time it is completed.
<i>Units</i>	None
<i>Notes</i>	None

2.35 Clockwise partial rotations of the animal's body

<i>Description</i>	Reports the number of times the animal's body rotated by at least the 'Partial rotation angle' defined in the protocol in a clockwise direction but didn't complete an entire rotation.
<i>Calculation method</i>	*** Feature still in beta ***
<i>Analysis across time</i>	This measure can be analysed across time. A partial rotation is deemed to occur at the time it is completed.
<i>Units</i>	None
<i>Notes</i>	None

2.36 Anti-clockwise partial rotations of the animal's body

<i>Description</i>	Reports the number of times the animal's body rotated by at least the 'Partial rotation angle' defined in the protocol in an anticlockwise direction but didn't complete an entire rotation.
<i>Calculation method</i>	*** Feature still in beta ***
<i>Analysis across time</i>	This measure can be analysed across time. A partial rotation is deemed to occur at the time it is completed.
<i>Units</i>	None
<i>Notes</i>	None

2.37 Absolute turn angle

<i>Description</i>	Reports the sum of the absolute angle between each movement vector of the animal.
<i>Calculation method</i>	A vector of movement from one position of the animal's centre point to the next is created. For each vector, the angle between it and the previous vector is calculated with anti-clockwise movement being negative and clockwise movement being positive (i.e. the angle is from -180° to 180°). The absolute value of this angle is summed for all the positions of the animal throughout the test or time period.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions within the specific time period.
<i>Units</i>	Degrees
<i>Notes</i>	From this measure, it is easy to use calculations to derive measures such as <i>Meander</i> and <i>Angular velocity</i> . The former is the <i>Absolute turn angle</i> divided by the <i>Total distance travelled</i> and the latter is the <i>Absolute turn angle</i> divided by the <i>Test duration</i> .

2.38 Absolute head turn angle

<i>Description</i>	Reports the cumulative absolute angle through which the animal's head moved. For example, if the animal moved its head 30° to the left and then moved its head 45° to the right, the absolute head turn angle would be 75° .
<i>Calculation method</i>	For each position of the animal's head a vector is created from the animal's centre point to the head. The angle between this vector and the same vector for the previous position of the animal's head is calculated and the absolute value of this angle is summed throughout the test.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	None

2.39 Clockwise head turn angle

<i>Description</i>	Reports the cumulative clockwise angle through which the animal's head moved. For example, if the animal moved its head 30° to the left and then moved its head 45° to the right, the clockwise head turn angle would be 45°.
<i>Calculation method</i>	For each position of the animal's head a vector is created from the animal's centre point to the head. The angle between this vector and the same vector for the previous position of the animal's head is calculate and if the angle is positive it is added to the sum of clockwise head turn angles for the test.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	None

2.40 Anti-clockwise head turn angle

<i>Description</i>	Reports the cumulative anti-clockwise angle through which the animal's head moved. For example, if the animal moved its head 30° to the left and then moved its head 45° to the right, the anti-clockwise head turn angle would be 30°.
<i>Calculation method</i>	For each position of the animal's head a vector is created from the animal's centre point to the head. The angle between this vector and the same vector for the previous position of the animal's head is calculate and if the angle is negative it is added to the sum of clockwise head turn angles for the test.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	None

2.41 Number of rears

<i>Description</i>	Reports the number of times the animal reared.
<i>Calculation method</i>	Counts the number of times the animal started to rear.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.42 Total time rearing

<i>Description</i>	Reports the total amount of time for which the animal was rearing.
--------------------	--

<i>Calculation method</i>	Sums the duration of each bout of rearing that occurred during the test. If the animal was rearing at the end of the test, then the last bout of rearing ends with the test end.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.43 Latency to first rear

<i>Description</i>	Reports the latency to the first time that the animal reared.
<i>Calculation method</i>	The time when the first bout of rearing started.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.44 Average duration of a rear

<i>Description</i>	Reports the average duration of the rearing bouts.
<i>Calculation method</i>	The result of Total time rearing divided by Number of rears.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.45 Maximum duration of a rear

<i>Description</i>	Reports the duration of the longest bout of rearing.
<i>Calculation method</i>	The duration of each bout of rearing is calculated and the longest bout is found.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.46 Minimum duration of a rear

<i>Description</i>	Reports the duration of the shortest bout of rearing.
<i>Calculation method</i>	The duration of each bout of rearing is calculated and the shortest bout is found.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

2.47 Path efficiency

<i>Description</i>	This measure represents an index of the efficiency of the path taken by the animal to get from the first position in the test to the last position. A value of 1 indicates perfect efficiency - the animal moved in a straight line - values less than 1 indicate decreasing efficiency.
<i>Calculation method</i>	The straight-line distance between the first position in the test and the last position is divided by the total distance travelled by the animal during the test.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the straight-line distance between the animal's position at the start of the time period or time segment and the animal's position at the end of the period or segment, divided by the distance the animal travelled during the period or segment. (Also see notes below.)
<i>Units</i>	None
<i>Notes</i>	<p>This measure is intended for use in water-maze experiments, but is available in all tests.</p> <p>This measure cannot be calculated if the animal passed through a hidden zone during the test, as the distance it will have travelled while hidden is not known. In such cases the measure's result will be shown as #N/A.</p> <p>Analysing this measure across time is usually not very meaningful, but can be helpful where a time marker is set at some specific moment in a test and the time marker is subsequently used to define a time period. In this situation the result for the time period can be used to determine path efficiency up to the moment the time marker was set.</p>

2.48 Number of line crossings

<i>Description</i>	Reports the number of times the animal's centre point moved from one area of the apparatus map to another - i.e. crossed the lines which constitute the map.
<i>Calculation method</i>	The apparatus is divided into unique areas by the apparatus map. For each animal position recorded in the experiment, the area which contains the animal's centre point is found. Each time this changes the measure's value is increased by 1.
<i>Analysis across time</i>	This measure can be analysed across time.

<i>Units</i>	None
<i>Notes</i>	This measure DOES NOT count transitions between ZONES. It counts transitions between areas of the apparatus map, irrespective of whether the areas are part of a zone. This measure is primarily intended to provide an easy way to measure 'grid line crossings' in a similar way to that commonly used manually - viz.: A regular grid is drawn on the apparatus and the experimenter counts the number of times the animal moves from one grid square to another. It's important to understand that ANY-maze uses the animal's centre point when calculating the measure and therefore it can be prone to 'spurious entries' if an animal straddles a line between two areas (i.e. by moving a very small amount, the animal can apparently cross a line many times). This problem can be overcome by setting one zone for each area, using the percentage of the animal that's in the zone to score zone entries, and then using a calculation to sum all the entries into these zones.

2.49 On/off inputs positive reversal

<i>Description</i>	The number of times the sequence in which on/off inputs were being activated changed from a decreasing sequence to an increasing one.
<i>Calculation method</i>	This measure is only applied to on/off inputs which have index values defined. As an input is activated, the system checks to see whether the index value of the newly activated input is greater or less than the index value of the previously active input. If the newly activated input has a lower index value, then the inputs are being activated in a decreasing sequence; whereas if it is higher, then they are being activated in an increasing sequence. Once the 'direction' of the sequence has been determined, then any change in direction is deemed to be a <i>reversal</i> . This measure counts the number of times the sequence changes from decreasing to increasing.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the number of positive reversals which occurred during the time period.
<i>Units</i>	None
<i>Notes</i>	This measure will only be available if the protocol includes two or more On/off inputs which have index values defined.

2.50 On/off inputs negative reversal

<i>Description</i>	The number of times the sequence in which on/off inputs were being activated changed from an increasing sequence to a decreasing one.
<i>Calculation method</i>	This measure is only applied to on/off inputs which have index values defined. As an input is activated, the system checks to see whether the index value of the newly activated input is greater or less than the index value of the previously active input. If the newly activated input has a lower index value, then the inputs are being activated in a decreasing sequence; whereas if it is higher, then they are being activated in an increasing sequence. Once the 'direction' of the sequence has been determined, then any change in direction is deemed to be a <i>reversal</i> . This measure counts the number of times the sequence changes from increasing to decreasing.

<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the number of negative reversals which occurred during the time period.
<i>Units</i>	None
<i>Notes</i>	This measure will only be available if the protocol includes two or more On/off inputs which have index values defined.

2.51 Tracking quality

<i>Description</i>	Classifies the tracking in a test as either <i>Good</i> or <i>Poor</i> .
<i>Calculation method</i>	The determination of tracking quality is based on two things: how many frames captured by the camera were successfully tracked, and the total amount of time that the animal was not tracked for. The actual calculation is complex as it takes into consideration the duration of each non-tracked period, with long periods being seen as much worse than short periods. For example, in a 2-minute test the animal might not be tracked for half a second twenty times during the test - however, ANY-maze would never have 'lost' the animal for more than half a second so the tracking errors this would cause would be minor. On the other hand if in a 2-minute test the animal was not tracked for a single period of 10 seconds then the <i>total</i> time the animal was 'lost' would be the same as in the first example, but in this case the tracking error could be quite significant (the animal could have moved quite a lot in 10 seconds). The calculation also takes into consideration the duration of the test. For example. if an animal is 'lost' for 30 seconds in a 2-minute test then that is quite significant, but if it was lost for 30 seconds in a 24 hour test then probably this would make no real difference to the test results.
<i>Units</i>	None
<i>Notes</i>	<p>ANY-maze's assessment of the tracking as <i>Good</i> or <i>Poor</i> is intended as a guide to the experimenter. You may, for example, wish to review tests ANY-maze classifies as having poor tracking.</p> <p>This measure can be useful when the option to remove jumps is switched on, as the removal of jumps can leave a test with long untracked periods, which will cause the test's tracking to be classified as <i>Poor</i>.</p>

2.52 Number of centre positions recorded

<i>Description</i>	Reports the number of positions of the centre of the animal that were tracked and recorded in the results. Note that ANY-maze may have <i>tracked</i> the centre of the animal in more positions than are reported here, because not all position are necessarily stored in the test's results - this depends on the setting for the maximum number of positions to record each second.
<i>Calculation method</i>	Counts the number of centre positions in the results.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the number of positions of the centre of the animal that fall within the time period.
<i>Units</i>	None
<i>Notes</i>	None

2.53 Number of head positions recorded

<i>Description</i>	Reports the number of positions of the head of the animal that were tracked and recorded in the results. Note that ANY-maze may have <i>tracked</i> the head of the animal in more positions than are reported here, because not all positions are necessarily stored in the test's results - this depends on the setting for the maximum number of positions to record each second.
<i>Calculation method</i>	Counts the number of head positions in the results.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the number of positions of the head of the animal that fall within the time period.
<i>Units</i>	None
<i>Notes</i>	The number of tail positions tracked is identical to the number of head positions tracked.

2.54 RAPC - Type 1 errors

<i>Description</i>	Reports the total number of <i>Type 1 errors</i> in the RAPC apparatus. A <i>Type 1 error</i> occurs when the animal tries to open a door that is latched shut.
<i>Calculation method</i>	The number of door 'openings' for all the doors in the RAPC apparatus is analysed (note that doors which are latched shut will still be registered as 'opening' when the animal pushes against the door, because the door will move a few millimetres). The last door opened in each chamber is, necessarily, the non-latched door; therefore the other doors must be latched. The number of openings of the latched doors is summed and this is the total number of <i>Type 1 errors</i> .
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the protocol includes 12 'switch inputs' with 'index' values of 1 - 12.

2.55 RAPC - Type 2 errors

<i>Description</i>	Reports the total number of <i>Type 2 errors</i> in the RAPC apparatus. A <i>Type 2 error</i> occurs when the animal opens a non-latched door but does not go through it into the next chamber.
<i>Calculation method</i>	The number of door 'openings' for all the doors in the RAPC apparatus is analysed (note that doors which are latched shut will still be registered as 'opening' when the animal pushes against the door, because the door will move a few millimetres). The last door opened in each chamber is, necessarily, the non-latched door; therefore the other doors must be latched. The number of openings of the non-latched door for each chamber less 1 is the number of <i>Type 2 errors</i> for that chamber. The sum for all the chambers is the total number of <i>Type 2 errors</i> .
<i>Units</i>	None

Notes This measure is only available if the protocol includes 12 'switch inputs' with 'index' values of 1 - 12.

2.56 RAPC - Door sequence

Description Reports the sequence of non-latched doors in the RAPC apparatus, where the doors in each chamber are numbered 1 through 3. Thus a value of 1321 would mean that door 1 between the first and second chamber was not latched, door 3 between the second and third chamber was not latched, and so on.

Calculation method The number of door 'openings' for all the doors in the RAPC apparatus is analysed (note that doors which are latched shut will still be registered as 'opening' when the animal pushes against the door because the door will move a few millimetres). The last door opened in each chamber is, necessarily, the non-latched door.

Units None

Notes This measure is only available if the protocol includes 12 'switch inputs' with 'index' values of 1 - 12.

3 Zone measures

3.1 Number of entries to the zone

<i>Description</i>	Counts the number of times the animal entered the zone.
<i>Calculation method</i>	Depends on the method used to detect zone entries - see Choosing how ANY-maze should detect entries into a zone for more details.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	None

3.2 Number of exits from the zone

<i>Description</i>	Reports the number of times the animal exited from a zone.
<i>Calculation method</i>	Depends on the method used to detect zone entries - see Choosing how ANY-maze should detect entries into a zone for more details.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	None

3.3 Number of entries of the animal's head into the zone

<i>Description</i>	Counts the number of times the animal's head entered the zone.
<i>Calculation method</i>	Counts the number of times the animal's head position changed from being outside the zone to being inside it.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those head positions within the time period.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if Head tracking is turned on.

3.4 Time in the zone

<i>Description</i>	Reports the total amount of time the animal spent in the zone.
<i>Calculation method</i>	Calculated by summing the duration of each visit to the zone where a visit starts at the time of a zone entry and ends at the time of a zone exit.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the amount of the period that the animal spent in the zone. For example, if an animal entered a zone at time 45 seconds and exited it at time 80 seconds, then for the time period 30-60 seconds the result would be 15 seconds.
<i>Units</i>	Seconds
<i>Notes</i>	None

3.5 Time the animal's head was in the zone

<i>Description</i>	Reports the total amount of time that the animal's head was in the zone.
<i>Calculation method</i>	Calculated by summing the duration of each visit of the animal's head to the zone where a visit starts at the time the animal's head entered the zone and ends at the time the animal's head exited the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the amount of the period that the animal's head spent in the zone. For example, if an animal's head entered a zone at time 45 seconds and exited it at time 80 seconds, then for the time period 30-60 seconds, the result would be 15 seconds.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if Head tracking is turned on.

3.6 Was first zone entered

<i>Description</i>	Reports whether the zone was the first zone the animal entered in the test.
<i>Calculation method</i>	Calculated by detecting the first zone entry in the test. This is affected by the <i>Don't score any results in this zone until the first 'true' entry</i> option on the <i>Zone entry settings</i> page. See Choosing how ANY-maze should detect entries into a zone for more details.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	The result of this measure is either YES or NO, therefore when analysed it will be treated as a 2 level nominal value - see Statistical tests included in ANY-maze.

3.7 Number of bouts of investigation of the zone

<i>Description</i>	Reports the total number of times the animal investigated the zone.
<i>Calculation method</i>	<p>When the animal starts investigating the zone an investigation bout begins; when it stops investigating the zone the bout ends. The number of bouts of investigation is equal to the number of bouts that began.</p> <p>The animal is deemed to be investigating the zone when its head is outside the zone but within the distance specified in the protocol of the zone's border, and, if specified in the protocol, the animal is oriented towards the zone and/or the animal's centre is not in the zone.</p>
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of bouts that began during the period.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.8 Total time investigating the zone

<i>Description</i>	Reports the total time that the animal spent investigating the zone.
<i>Calculation method</i>	When the animal starts investigating the zone a note is made of the test clock; this is subtracted from the test clock when the animal stops investigating the zone - giving the duration of the investigation bout. The duration of all bouts is summed to give the total time investigating the zone. The animal is deemed to be investigating the zone when its head is outside the zone but within the distance specified in the protocol of the zone's border, and, if specified in the protocol, the animal is oriented towards the zone and/or the animal's centre is not in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is total time the animal was investigating the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.9 Was first zone investigated

<i>Description</i>	Reports whether the zone was the first zone the animal investigated during the test.
<i>Calculation method</i>	Calculated by detecting the first bout of investigation of a zone in the test.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the zone is an investigation zone. The result of this measure is either YES or NO, therefore when analysed it will be treated as a 2 level nominal value - see Statistical tests included in ANY-maze.

3.10 Longest bout of investigation of the zone

<i>Description</i>	Reports the duration of the longest bout of investigation of the zone.
<i>Calculation method</i>	The animal is deemed to be investigating the zone when its head is outside the zone but within the distance specified in the protocol of the zone's border, and, if specified in the protocol, the animal is oriented towards the zone and/or the animal's centre is not in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the duration of the longest bout of investigation of the zone in the period. If the animal is investigating the zone at the start of a time period then, for the purposes of analysis across time of this measure, the bout starts at the start of the period; likewise if the animal is investigating the zone at the end of the time period the bout ends at the end of the time period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.11 Shortest bout of investigation of the zone

<i>Description</i>	Reports the duration of the shortest bout of investigation of the zone.
<i>Calculation method</i>	The animal is deemed to be investigating the zone when its head is outside the zone but within the distance specified in the protocol of the zone's border, and, if specified in the protocol, the animal is oriented towards the zone and/or the animal's centre is not in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the duration of the shortest bout of investigation of the zone in the period. If the animal is investigating the zone at the start of a time period then, for the purposes of analysis across time of this measure, the bout starts at the start of the period; likewise if the animal is investigating the zone at the end of the time period the bout ends at the end of the time period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.12 Average duration of investigation of the zone

<i>Description</i>	Reports the average duration of bouts of investigation of the zone.
<i>Calculation method</i>	Calculated by dividing the Total time investigating the zone by the Number of bouts of investigation of the zone. If the animal did not investigate the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is Time spent investigating the zone during the time period divided by the Number of bouts of investigation of the zone in the time period. If the animal did not investigate the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone. This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

3.13 List of the duration of each investigation of the zone

<i>Description</i>	Reports a comma-separated list of the duration of each bout of investigation of the zone. For example, if the animal investigated the zone three times during the test, for 1 second on the first occasion and for 20 seconds on the second and third occasions then the list would be '1.0, 20.0, 20.0'.
<i>Calculation method</i>	The time that the animal starts investigating the zone is noted; then, when the animal stops investigating the zone, the duration of the bout of investigation is calculated and added to the list. If the animal is investigating the zone at the end of the test, the time from the start of the investigation to the end of the test is used as the duration of the last visit.
<i>Analysis across time</i>	This measure cannot be analysed across time.

<i>Units</i>	The duration of all bouts of investigation is reported in seconds.
<i>Notes</i>	<p>This measure is only available if the zone is an investigation zone.</p> <p>When included on the Data page, this measure will show all the bouts of investigation of the zone in a single cell. If the spreadsheet is saved in CSV format and then opened in (for example) Microsoft Excel, then the visits will be listed in individual cells.</p> <p>The length of the list is limited to 8192 characters, but usually at least 1,000 bouts of investigation will be listed before the limit is reached.</p>

3.14 Number of partial exits from the zone

<i>Description</i>	Reports the number of times the animal partially exited from the zone.
<i>Calculation method</i>	A partial exit occurs when more than the partial exit percentage (as specified in the Zone's entry and exit settings in the protocol) of the animal exits the hidden zone and the animal subsequently fully returns into the zone (if it fully exits the zone then a partial exit is NOT scored). This measure counts the number of times this occurs.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the number of partial exits that <i>started</i> during the period.
<i>Units</i>	None
<i>Notes</i>	<p>Only available for hidden zones that are set to score partial exits.</p> <p>Partial exits can only be scored in tests that were performed using ANY-maze version 7.20 or above. Provided the test was performed with a suitable version, then scoring of partial exits can be switched on after the tests have been performed, i.e. you don't need to have chosen to score partial exits at the time that the test was run.</p>

3.15 Time partially exited from the zone

<i>Description</i>	Reports the total amount of time the animal is partially exited from the zone.
<i>Calculation method</i>	A partial exit starts when more than the partial exit percentage (as specified in the Zone's entry and exit settings in the protocol) of the animal exits the hidden zone, and ends when the animal subsequently fully returns into the zone (if it fully exits the zone then a partial exit is NOT scored). This measure calculates the time between the start and the end of each partial exit and sums them to give the total amount of time partially exited from the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the amount of time the animal was partially exited from the zone during the period. If the animal is partially exited from the zone at the start of the period, then the result includes the time until the animal returns into the zone (or the period ends). This means that the count of partial exits from the zone during a time period can be zero, while the time partially exited from the zone during the same period is non-zero.
<i>Units</i>	Seconds

Notes

Only available for hidden zones that are set to score partial exits.
Partial exits can only be scored in tests that were performed using ANY-maze version 7.20 or above. Provided the test was **performed** with a suitable version, then scoring of partial exits can be switched on after the tests have been performed, i.e. you don't need to have chosen to score partial exits at the time that the test was run.

3.16 Distance travelled in the zone

<i>Description</i>	Reports the distance the animal travelled while in the zone.
<i>Calculation method</i>	Calculated by summing the distance travelled during each visit to the zone. A visit starts when an animal enters the zone and ends when it exits the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the distance travelled within the zone during that time period.
<i>Units</i>	Metres
<i>Notes</i>	As the position of the animal prior to a zone entry must, by definition, be outside the zone and the position after it enters the zone must be inside, the distance between the two positions will be partly outside and partly inside the zone. ANY-maze adds all this distance to the distance travelled in the zone the animal's <u>leaving</u> . Although this can lead to inaccuracies, they are generally very small because: a) ANY-maze detects many positions per second so the distance between any two positions is usually very small; b) Any small distance 'lost' when the animal enters a zone entry is usually counterbalanced by a small distance which is 'gained' when it leaves the zone - see figure 3.




Figure 3. When the zone change at point 'X' is detected, ANY-maze adds all the distance 'A' to the distance travelled in the 'Outer zone'. When the second zone change at point 'Y' is detected, ANY-maze adds all the distance 'B' to the distance travelled in the 'Inner zone'. [Note: Distances A and B have been exaggerated in this diagram to aid explanation.]

In some situations, tracks can have small oscillations in them which tend to generate unrepresentatively large values for distance travelled. This occurs most often when an animal travels slowly while moving its body a lot - for example, while exploring an open field. To overcome this, ANY-maze uses an adaptive smoothing algorithm to attenuate these oscillations when calculating distance travelled - see figure 4. Note: The definition of what's a small oscillation is based on the animal's size.




Figure 4. Measuring the length of the actual track (shown in orange) would yield an unrepresentatively large value for distance travelled. ANY-maze uses a 'smoothed' track (shown in green) to better estimate the true distance travelled. [Note: The oscillations in this track have been exaggerated to aid explanation.]

3.17 Distance travelled until first entry into the zone

<i>Description</i>	Reports the distance travelled by the animal up to its first entry into the specified zone.
<i>Calculation method</i>	The distance the animal travels is summed until it enters the zone. If the animal doesn't enter the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	This measure is affected by the <i>Don't score any results in this zone until the first 'true' entry</i> option on the <i>Zone entry settings</i> page. See <i>Choosing how ANY-maze should detect entries into a zone</i> for more details.

3.18 Distance travelled by the animal's head in the zone

<i>Description</i>	Reports the distance the animal's head travelled while the head was in the zone.
<i>Calculation method</i>	Calculated by summing the distance travelled by the animal's head during each visit to the zone. A visit starts when the animal's head enters the zone and ends when it exits the zone.

<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	As the position of the animal's head prior to a zone entry must, by definition, be outside the zone and the position after it enters the zone must be inside, the distance between the two positions will be partly outside and partly inside the zone. ANY-maze adds all this distance to the distance travelled by the animal's head in the zone that the animal is <i>leaving</i> . Although this can lead to inaccuracies, they are generally very small because: a) ANY-maze detects many positions per second so the distance between any two positions is usually very small; b) Any small distance 'lost' when the animal enters a zone entry is usually counterbalanced by a small distance which is 'gained' when it leaves the zone. This measure is only available if Head tracking is turned on.

3.19 Distance travelled while investigating the zone

<i>Description</i>	Reports the distance the animal travelled while investigating the zone.
<i>Calculation method</i>	Calculated by summing the distance travelled during each bout of investigation of the zone. A visit starts when an animal starts investigating the zone and ends when it stops investigating the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the distance travelled while investigating the zone during that time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.20 Distance travelled before first investigation of the zone

<i>Description</i>	Reports the distance travelled by the animal up to the point that it starts investigating the specified zone.
<i>Calculation method</i>	The distance the animal travels is summed until it starts investigating the zone. If the animal doesn't investigate the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.21 Latency to first entry to the zone

<i>Description</i>	Reports the length of time which elapsed before the animal entered the zone for the first time.
<i>Calculation method</i>	This value is the time at which the first zone entry occurred. If the animal doesn't enter the zone during the test, then the result is undefined.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the first entry. If the animal doesn't enter the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the <i>Don't score any results in this zone until the first 'true' entry</i> option on the <i>Zone entry settings</i> page. See Choosing how ANY-maze should detect entries into a zone for more details. This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.

3.22 Latency to first exit from the zone

<i>Description</i>	Reports the length of time which elapsed before the animal exited the zone for the first time.
<i>Calculation method</i>	This value is the time at which the first zone exit occurred. If the animal doesn't exit the zone during the test then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the first exit. If the animal doesn't exit the zone during the time period then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the <i>Don't score any results in this zone until the first 'true' entry</i> option on the <i>Zone entry settings</i> page. See Choosing how ANY-maze should detect entries into a zone for more details. This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.

3.23 Latency to last entry to the zone

<i>Description</i>	Reports the length of time that elapsed up to the moment when the animal made its last entry into the zone during the test.
<i>Calculation method</i>	This value is updated at the moment of each entry into the zone. The value at the end of the test is, necessarily, the latency to the last zone entry. If the animal doesn't enter the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the last entry during that period. If the animal doesn't enter the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds This measure is particularly useful in water-maze tests to report the time taken by the animal to find a platform zone. In such tests, it's common to only consider the animal as finding the platform if it remains on it for a certain period, for example 5 seconds. In this case, it's possible that the animal will enter the platform zone a number of times and therefore the time to 'find' the platform will be the latency to the <i>last</i> entry to the zone.

This measure is affected by the option to *Use the test duration as the latency for events which don't occur* in the Analysis element.

3.24 Latency to the first entry of the animal's head into the zone

<i>Description</i>	Reports the length of time which elapsed before the animal's head entered the zone for the first time.
<i>Calculation method</i>	This value is the time at which the animal's head first entered the zone. If the animal's head doesn't enter the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the first entry of the animal's head into the zone. If the animal's head doesn't enter the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if Head tracking is turned on. This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.

3.25 Latency to the first exit of the animal's head from the zone

<i>Description</i>	Reports the length of time which elapsed before the animal's head exited the zone for the first time.
<i>Calculation method</i>	This value is the time at which the animal's head first exited the zone. If the animal doesn't exit the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the first exit of the animal's head from the zone. If the animal's head doesn't exit the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if Head tracking is turned on. This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.

3.26 Latency to first investigation of the zone

<i>Description</i>	Reports the time that elapsed in the test before the animal investigated the zone for the first time.
<i>Calculation method</i>	The value of the test clock when the first bout of investigation began. The animal is deemed to be investigating the zone when its head is outside the zone but within the distance specified in the protocol of the zone's border, and, if specified in the protocol, the animal is oriented towards the zone and/or the animal's centre is not in the zone.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the time of the first bout of investigation in the period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.27 Latency to end of first investigation of the zone

<i>Description</i>	Reports the length of time which elapsed before the animal's first investigation of the zone ended.
<i>Calculation method</i>	This value is the time at which the animal ended its first bout of investigation of the zone. If the animal doesn't investigate the zone during the test then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time from the start of the period to the end of the animal's first bout of investigation of the zone. If the animal doesn't investigate the zone during the time period then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone. This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.

3.28 Average speed in the zone

<i>Description</i>	Reports the average speed of the animal while it was in the zone.
<i>Calculation method</i>	Calculated by dividing the Distance travelled in the zone by the Time in the zone. If the animal was never in the zone during the test then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Distance travelled in the zone during the time period divided by the Time spent in the zone during the time period. If the animal was never in the zone during the time period then the result is undefined.
<i>Units</i>	Metres per second
<i>Notes</i>	If you want to know the average speed in the zone while mobile (i.e. ignoring periods when the animal was stationary), then use a calculation of <i>Distance travelled in the zone / Time mobile in the zone</i> . This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

3.29 Maximum speed in the zone

<i>Description</i>	Reports the maximum speed of the animal while in the zone.
--------------------	--

<i>Calculation method</i>	The speed of the animal between positions within the zone is calculated and the maximum speed is found. If the animal doesn't enter the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. If the animal doesn't enter the zone during the time period, then the result is undefined.
<i>Units</i>	Metres per second
<i>Notes</i>	The calculation of maximum speed does not use <i>successive</i> positions but instead requires that the animal move at least a minimum distance (which is based on the animal's size) and the speed to cover <i>this</i> distance is calculated. This method of calculation is used to avoid reporting the speed of movements that don't constitute locomotion of the animal. For example, if an animal scratches, its centre point may oscillate rapidly but this will not be reported as the animal's maximum speed.

3.30 Average speed while investigating the zone

<i>Description</i>	Reports the average speed of the animal while it was investigating the zone.
<i>Calculation method</i>	Calculated by dividing the Distance travelled while investigating the zone by the Time spent investigating the zone. If the animal never investigated the zone during the test then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Distance travelled while investigating the zone during the time period divided by the Time spent investigating the zone during the time period. If the animal did not investigate the zone during the time period then the result is undefined.
<i>Units</i>	Metres per second
<i>Notes</i>	This measure is only available if the zone is an investigation zone. This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

3.31 Longest visit to the zone

<i>Description</i>	Reports the duration of the longest single visit to the zone.
<i>Calculation method</i>	Each zone visit is delimited by a zone entry and a zone exit - the time between the two is the duration of the visit. The duration of each visit is calculated and the largest value is found. If the animal was never in the zone during the test then the result is zero.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the duration of the longest visit to the zone during the time period. If the animal spent the entire period in the zone, then the result will be the duration of the time period itself.
<i>Units</i>	Seconds
<i>Notes</i>	None

3.32 Shortest visit to the zone

<i>Description</i>	Reports the duration of the shortest single visit to the zone.
<i>Calculation method</i>	Each zone visit is delimited by a zone entry and a zone exit - the time between the two is the duration of the visit. The duration of each visit is calculated and the smallest value is found. If the animal was never in the zone during the test, then the result is zero.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the duration of the shortest visit to the zone during the time period.
<i>Units</i>	Seconds
<i>Notes</i>	None

3.33 Average duration of visit to the zone

<i>Description</i>	Reports the average duration of visits to the zone.
<i>Calculation method</i>	Calculated by dividing the Time spent in the zone by the Number of entries to the zone. If the animal was never in the zone during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is Time spent in the zone during the time period divided by the Number of entries to the zone in the time period. If the animal was never in the zone during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

3.34 List of the duration of each visit to the zone

<i>Description</i>	Reports a comma-separated list of the duration of each visit to the zone. For example, if the animal visited the zone three times during the test, for 1 second on the first occasion and for 20 seconds on the second and third occasions then the list would be '1.0, 20.0, 20.0'.
<i>Calculation method</i>	The time of the animal's entry to the zone is noted; then, when the animal exits the zone, the duration of the visit is calculated and added to the list. If the animal is in the zone at the end of the test, the time from the entry to the end of the test is used as the duration of the last visit.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	The duration of all visits is reported in seconds.
<i>Notes</i>	When included on the Data page, this measure will show all the visits to the zone in a single cell. If the spreadsheet is saved in CSV format and then opened in (for example) Microsoft Excel, then the visits will be listed in individual cells.

The length of the list is limited to 8192 characters, but usually at least 1,000 visits will be listed before the limit is reached.

3.35 Time mobile in the zone

<i>Description</i>	Reports the total time that the animal was mobile in the zone.
<i>Calculation method</i>	Calculated by subtracting the Time immobile in the zone from the Time spent in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Time spent in the zone during the period minus the Time immobile in the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

3.36 Time immobile in the zone

<i>Description</i>	Reports the total time that the animal was immobile in the zone.
<i>Calculation method</i>	Sums the duration of each immobile episode in the zone - see notes for the definition of an immobile episode.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the sum of the duration of each immobile episode in the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	<p>An immobile episode in the zone starts when the animal becomes immobile after being mobile or when the animal is already immobile and enters the zone. Generally, it's unlikely that an immobile animal will enter a zone, because to enter the zone it will probably have to be mobile. Nevertheless, if the animal is immobile right on the border of a zone, it could enter the zone by moving very slightly but not by enough to end the immobile episode.</p> <p>An immobile episode in the zone ends when the animal becomes mobile or when it leaves the zone.</p> <p>The definition of immobility depends on the protocol - see Immobility detection.</p>

3.37 Immobile episodes in the zone

<i>Description</i>	Reports the number of times the animal became immobile while in the zone.
<i>Calculation method</i>	Counts the number of times the animal changed from being mobile to being immobile while in the zone. If an <i>immobile</i> animal enters a zone (see note below), then the entry will be considered to start a new immobile episode in the zone, i.e. the count of immobile episodes in the zone will be incremented. This means that the sum of the immobile episodes in all the zones may be greater than the number of immobile episodes in the apparatus as a whole.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of immobile episodes that started in the zone during the time period. If an animal is already immobile in the zone at the start of the time period, then a new immobile episode is NOT counted for the period. This means that it's possible to have a period for which the result of this measure is zero but the result for the Time immobile in the zone is not zero.
<i>Units</i>	None
<i>Notes</i>	An immobile episode in the zone starts when the animal becomes immobile after being mobile or when the animal is already immobile and enters the zone. Generally, it's unlikely that an immobile animal will enter a zone, because to enter the zone it will probably have to be mobile. Nevertheless if the animal is immobile right on the border of a zone it could enter the zone by moving very slightly but not by enough to end the immobile episode. An immobile episode in the zone ends when the animal becomes mobile or when it leaves the zone. The definition of immobility depends on the protocol - see Immobility detection.

3.38 Time mobile while investigating the zone

<i>Description</i>	Reports the total time that the animal was mobile while investigating the zone.
<i>Calculation method</i>	Calculated by subtracting the Time immobile while investigating the zone from the Time spent investigating the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Time spent investigating the zone during the period minus the Time immobile while investigating the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone. The definition of immobility depends on the protocol - see Immobility detection.

3.39 Time immobile while investigating the zone

<i>Description</i>	Reports the total time that the animal was immobile while investigating the zone.
<i>Calculation method</i>	Sums the duration of each immobile episode while investigating the zone - see notes for the definition of an immobile episode.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the sum of the duration of each immobile episode while the animal was investigating the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone. An immobile episode while investigating the zone starts when the animal becomes immobile after being mobile or when the animal is already immobile and starts investigating the zone. Generally, it's unlikely that an immobile animal will start investigating a zone, because to start investigating the zone it will probably have to

move. Nevertheless, if the animal is immobile right on the border of the zone's investigation area, it could start investigating the zone by moving very slightly but not by enough to end the immobile episode.

An immobile episode while investigating the zone ends when the animal becomes mobile or when it stops investigating the zone.

The definition of immobility depends on the protocol - see Immobility detection.

3.40 Immobile episodes while investigating the zone

<i>Description</i>	Reports the number of times the animal became immobile while investigating the zone.
<i>Calculation method</i>	Counts the number of times the animal changed from being mobile to being immobile while investigating the zone. If an <i>immobile</i> animal starts investigating a zone (see note below), then the entry will be considered to start a new immobile episode while investigating the zone, i.e. the count of immobile episodes while investigating the zone will be incremented. This means that the sum of the immobile episodes while investigating all the zones may be greater than the number of immobile episodes in the apparatus as a whole.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of immobile episodes that started while the animal was investigating the zone during the time period. If an animal is already immobile while investigating the zone at the start of the time period, then a new immobile episode is NOT counted for the period. This means that it's possible to have a period for which the result of this measure is zero but the result for the Time immobile while investigating the zone is not zero.
<i>Units</i>	None
<i>Notes</i>	<p>This measure is only available if the zone is an investigation zone.</p> <p>An immobile episode while investigating the zone starts when the animal becomes immobile after being mobile or when the animal is already immobile and starts investigating the zone. Generally, it's unlikely that an immobile animal will start investigating a zone, because to enter the zone it will probably have to move. Nevertheless if the animal is immobile right on the border of the zone's investigation area, it could start investigating the zone by moving very slightly but not by enough to end the immobile episode.</p> <p>An immobile episode while investigating the zone ends when the animal becomes mobile or when it stops investigating the zone.</p> <p>The definition of immobility depends on the protocol - see Immobility detection.</p>

3.41 Time active in the zone

<i>Description</i>	Reports the total time that the animal was active in the zone.
<i>Calculation method</i>	Calculated by subtracting the Time inactive in the zone from the Time in the zone.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Time spent in the zone during the period minus the Time inactive in the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	An animal is active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - grooming for example.

3.42 Time inactive in the zone

<i>Description</i>	Reports the total time that the animal was inactive in the zone.
<i>Calculation method</i>	Sums the duration of each inactive episode in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the sum of the duration of each inactive episode in the zone during the period.
<i>Units</i>	Seconds
<i>Notes</i>	<p>An inactive episode in the zone starts when the animal becomes inactive after being active or when the animal enters the zone and is already inactive. An inactive episode in the zone ends when an animal becomes active or when it leaves the zone.</p> <p>Inactivity is defined as NOT activity. An animal is defined to be active if it is either mobile OR it's performing some other behaviour which has been specified as an activity - for example, grooming. If the protocol specifies that immobility should not be detected, then activity analysis will be based purely on the performance of other behaviours.</p>

3.43 Inactive episodes in the zone

<i>Description</i>	Reports the total number of times the animal became inactive while in the zone.
<i>Calculation method</i>	Counts the number of times the animal changed from being active to being inactive while in the zone. If an <i>inactive</i> animal enters a zone, then the entry will be considered to start a new inactivity episode in the zone, i.e. the count of inactive episodes in the zone will be incremented. This means that the sum of the inactive episodes in all the zones may be greater than the number of inactive episodes in the apparatus as a whole.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of inactive episodes that started in the zone during the time period. If an animal is already inactive in the zone at the start of the time period, then a new inactive episode is NOT counted for the period. This means that it's possible to have a period for which the result of this measure is zero but the result for the Time inactive in the zone is not zero.
<i>Units</i>	None
<i>Notes</i>	An inactive episode in the zone starts when the animal becomes inactive after being active or when the animal enters the zone and is already inactive. An inactive episode in the zone ends when an animal becomes active or when it leaves the zone.

3.44 Initial distance from the zone

<i>Description</i>	Reports the distance from the animal's first position in the test to the zone
<i>Calculation method</i>	The straight line distance from the first position of the animal (see notes) to the nearest point of the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The value is the distance of the animal from the zone at the start of the time period.
<i>Units</i>	Metres
<i>Notes</i>	The position of the animal is either the animal's centre or the point of the animal that is closest to the zone. Which of these is used depends on whether or not the zone entry settings use the animal area to determine zone entries; when they do, the distance from the animal to the zone is based on the distance from the point of the animal that is closest to the zone; otherwise, the distance is based on the centre of the animal.

3.45 Average distance from the zone

<i>Description</i>	Reports the average distance from the animal to the zone when the animal is <i>outside</i> the zone.
<i>Calculation method</i>	<p>ANY-maze calculates the distance from the animal to the zone for every position of the animal that is outside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand if zone entries are based on the centre of the animal, then the distance to the zone will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest part of the zone.</p> <p>Having determined the distance from the zone, ANY-maze maintains a cumulative sum of each distance multiplied by the time the animal remained at that distance. The final result for the average distance from the zone is this cumulative sum divided by the total duration of the test or period.</p> <p>The reason the system works this way is best explained using an example. Imagine the animal was 50cm from a zone and remained there for 55 seconds; it then moved to be 30cm from the zone and remained there for 5 seconds; the test then ended. Just taking the average of the two distances would imply that the average distance from the zone was 40cm, but this is very misleading as the animal spent almost the entire test 50cm from the zone. Instead, ANY-maze would calculate the average distance as $[(50 \times 55) + (30 \times 5)] / 60 = 48.33\text{cm}$. Effectively, the system weights the distances depending on how long the animal remained there.</p>

<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal spends the entire duration of the test (or of a time period) inside the zone, then the result will be zero (i.e. the animal was no distance from the zone).

3.46 Maximum distance from the zone

<i>Description</i>	Reports the maximum distance from the animal to the zone when the animal is <i>outside</i> the zone.
<i>Calculation method</i>	ANY-maze calculates the distance from the animal to the zone for every position of the animal that is outside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal, or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest part of the zone. Having determined the distance from the animal to the zone, the system simply notes the maximum value during the test or time period.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the maximum distance considering just those positions of the animal that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal spends the entire duration of the test (or of a time period) inside the zone, then the result will be zero (i.e. the animal was no distance from the zone).

3.47 Minimum distance from the zone

<i>Description</i>	Reports the minimum distance from the animal to the zone when the animal is <i>outside</i> the zone.
<i>Calculation method</i>	ANY-maze calculates the distance from the animal to the zone for every position of the animal that is outside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal, or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest part of the zone. Having determined the distance from

the animal to the zone, the system simply notes the minimum value during the test or time period.

Analysis across time This measure can be analysed across time. The result is the minimum distance considering just those positions of the animal that fall within the time period.

Units Metres

Notes If the animal enters the zone, then this value is automatically set to zero.

3.48 Cumulative distance from the zone

Description Reports the sum of the product of the distance from the zone and the time at that distance.

Calculation method For every position of the animal, this calculates the distance from the zone multiplied by the time the animal stayed at that position. The final result is the sum of all these values.

Analysis across time This measure can be analysed across time.

Units Metres · seconds

Notes This value represents the area under the curve of a graph of distance from zone vs. time.

3.49 Average distance of the animal's head from the zone

Description Reports the average distance from the animal's head to the zone when the animal is *outside* the zone.

Calculation method ANY-maze calculates the distance from the animal's head to the closest point on the zone border for every position of the animal's head that is outside the zone. The system maintains a cumulative sum of each distance multiplied by the time the animal remained at that distance. The final result for the average distance of the animal's head from the zone is this cumulative sum divided by the total duration of the test or period.

The reason the system works this way is best explained using an example. Imagine the animal's head was 50cm from a zone and remained there for 55 seconds; it then moved to be 30cm from the zone and remained there for 5 seconds; the test then ended. Just taking the average of the two distances would imply that the average distance from the zone was 40cm, but this is very misleading as the animal spent almost the entire test 50cm from the zone. Instead, ANY-maze would calculate the average distance as $[(50 \times 55) + (30 \times 5)] / 60 = 48.33\text{cm}$. Effectively, the system weights the distances depending on how long the animal remained there.

Analysis across time This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.

Units Metres

Notes This measure is only available if Head tracking is turned on.

If the animal's head spends the entire duration of the test (or of a time period) inside the zone, then the result will be zero (i.e. the animal's head was no distance from the zone).

3.50 Maximum distance of the animal's head from the zone

<i>Description</i>	Reports the maximum distance from the animal's head to the zone when the animal is <i>outside</i> the zone.
<i>Calculation method</i>	For each position of the animal's head that is outside the zone, ANY-maze calculates the distance from the head to the closest point on the zone border. The maximum such distance is found.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if Head tracking is turned on. If the animal's head spends the entire duration of the test (or of a time period) inside the zone, then the result will be zero (i.e. the animal's head was no distance from the zone).

3.51 Minimum distance of the animal's head from the zone

<i>Description</i>	Reports the minimum distance from the animal's head to the zone when the animal is <i>outside</i> the zone.
<i>Calculation method</i>	For each position of the animal's head that is outside the zone, ANY-maze calculates the distance from the head to the closest point on the zone border. The minimum such distance is found.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if Head tracking is turned on. If the animal's head enters the zone, then this value will be zero.

3.52 Average distance to the zone border

<i>Description</i>	Reports the average distance from the animal to the border of the zone when the animal is <i>inside</i> the zone.
<i>Calculation method</i>	ANY-maze calculates the distance from the animal to the zone border for every position of the animal that is inside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal then the calculation of the distance from the animal to the zone border will also be based

on the entire area of the animal - specifically the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand if zone entries are based on the centre of the animal then the distance to the zone border will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest border of the zone.

Having determined the distance to the zone border ANY-maze maintains a cumulative sum of each distance multiplied by the time the animal remained at that distance. The final result for the average distance to the zone border is this cumulative sum divided by the total duration of the test or period.

The reason the system works this way is best explained using an example. Imagine the animal was 20cm from a zone border and remained there for 55 seconds; it then moved to be 10cm from the zone border and remained there for 5 seconds; the test then ended. Just taking the average of the two distances would imply that the average distance from the zone border was 15cm, but this is very misleading as the animal spent almost the entire test 20cm from the border. Instead, ANY-maze would calculate the average distance as $[(20 \times 55) + (10 \times 5)] / 60 = 19.16\text{cm}$. Effectively, the system weights the distances depending on how long the animal remained there.

<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal never enters the zone, then the result is either undefined or zero depending on the setting <i>Use zero as the result for undefined averages</i> in the Analysis options element.

3.53 Maximum distance to the zone border

<i>Description</i>	Reports the maximum distance from the animal to the border of the zone when the animal is <i>inside</i> the zone.
<i>Calculation method</i>	ANY-maze calculates the distance from the animal to the zone border for every position of the animal that is inside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone border will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone border will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest border of the zone. Having determined the distance to the border, the system simply notes the maximum value during the test or time period.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the maximum distance considering just those positions of the animal that fall within the time period.
<i>Units</i>	Metres

Notes If the animal never enters the zone, then the result is undefined.

3.54 Minimum distance to the zone border

Description Reports the minimum distance from the animal to the border of the zone when the animal is *inside* the zone.

Calculation method ANY-maze calculates the distance from the animal to the zone border for every position of the animal that is inside the zone. Exactly how this is done depends on whether zone entries (*sic*) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone border will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone border will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest border of the zone. Having determined the distance to the border, the system simply notes the minimum value during the test or time period.

Analysis across time This measure can be analysed across time. The result is the minimum distance considering just those positions of the animal that fall within the time period.

Units Metres

Notes If the animal never enters the zone, then the result is undefined. If the animal exits the zone, this value is automatically set to zero.

3.55 Average distance from the animal's head to the zone border

Description Reports the average distance from the animal's head to the border of the zone when the animal is *inside* the zone.

Calculation method ANY-maze calculates the distance from the animal's head to the closest point on the zone border for every position of the animal's head that is inside the zone. The system maintains a cumulative sum of each distance multiplied by the time the animal remained at that distance. The final result for the average distance to the zone border is this cumulative sum divided by the total duration of the test or period.

The reason the system works this way is best explained using an example. Imagine the animal's head was 20cm from a zone border and remained there for 55 seconds; it then moved to be 10cm from the zone border and remained there for 5 seconds; the test then ended. Just taking the average of the two distances would imply that the average distance from the animal's head to the zone border was 15cm, but this is very misleading as the animal's head spent almost the entire test 20cm from the border. Instead, ANY-maze would calculate the average distance as $[(20 \times 55) + (10 \times 5)] / 60 = 19.16\text{cm}$. Effectively, the system weights the distances depending on how long the animal remained there.

<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal's head never enters the zone, then the result is either undefined or zero depending on the setting <i>Use zero as the result for undefined averages</i> in the Analysis options element. This measure is only available if Head tracking is turned on.

3.56 Maximum distance from the animal's head to the zone border

<i>Description</i>	Reports the maximum distance from the animal's head to the border of the zone when the animal is <i>inside</i> the zone.
<i>Calculation method</i>	For each position of the animal's head that is inside the zone, ANY-maze calculates the distance from the head to the closest point on the zone border. The maximum such distance is found.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal's head never enters the zone, then the result is undefined. This measure is only available if Head tracking is turned on.

3.57 Minimum distance from the animal's head to the zone border

<i>Description</i>	Reports the minimum distance from the animal's head to the border of the zone when the animal is <i>inside</i> the zone.
<i>Calculation method</i>	For each position of the animal's head that is inside the zone, ANY-maze calculates the distance from the head to the closest point on the zone border. The minimum such distance is found.
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal's head that fall within the time period.
<i>Units</i>	Metres
<i>Notes</i>	If the animal's head never enters the zone, then the result is undefined. This measure is only available if Head tracking is turned on.

3.58 Time getting closer to the zone

<i>Description</i>	Reports the total amount of time that the animal was outside the zone and was getting closer to it.
<i>Calculation method</i>	ANY-maze calculates the distance from the animal to the zone for every position of the animal that is outside the zone. Exactly how this is done depends on whether zone entries (<i>sic</i>) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If

entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest part of the zone.

Having calculated the distance to the zone, ANY-maze compares it to the previous distance to the zone; if it is less, then the animal is getting closer to the zone and the time from the previous position to this one is added to the total time getting closer to the zone. Note that very small movements of the animal will be ignored by the system (the definition of very small being based on the animal's size).

Analysis across time This measure can be analysed across time. The result is based on just those positions of the animal that fall within the time period.

Units Seconds

Notes This measure seems very similar to the Time moving towards the zone, but is calculated quite differently (see the definition of Time moving towards the zone for details on how it's calculated). The principal difference is that this measure relates to the animal's distance from the zone, whereas the *Time moving towards the zone* measure relates to the animal's heading. A good example of this difference is when a zone is in the form of a ring. In this case, if the animal is moving inside the ring it would always be 'moving towards' the zone as the zone surrounds it, but it could still either be getting closer or further from the zone. (In fact, it would always be moving both towards and away from the zone, so because of this ambiguity, ANY-maze simply wouldn't score either of these measures.)

3.59 Time getting further away from the zone

Description Reports the total amount of time that the animal was outside the zone and was getting further away from it.

Calculation method ANY-maze calculates the distance from the animal to the zone for every position of the animal that is outside the zone. Exactly how this is done depends on whether zone entries (*sic*) are set to use the entire area of the animal or the animal's centre point (see Choosing how ANY-maze should detect entries into a zone for details). If entries are based on the entire area of the animal, then the calculation of the distance from the animal to the zone will also be based on the entire area of the animal - specifically, the system will use the distance from the point on the animal's edge that is closest to the zone border; on the other hand, if zone entries are based on the centre of the animal then the distance to the zone will also be based on the centre of the animal - i.e. the distance to the zone will simply be the distance from the centre to the nearest part of the zone.

Having calculated the distance to the zone, ANY-maze compares it to previous distance to the zone; if it is greater, then the animal is getting further away from the zone and the time from the previous position to this one is added to the total time getting further away from the zone. Note that very small movements of the

animal will be ignored by the system (the definition of very small being based on the animal's size).

<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions of the animal that fall within the time period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure seems very similar to the Time moving away from the zone, but is calculated quite differently (see the definition of Time moving away from the zone for details on how it's calculated). The principal difference is that this measure relates to the animal's distance from the zone, whereas the <i>Time moving away from the zone</i> measure relates to the animal's heading. A good example of this difference is when a zone is in the form of a ring. In this case, if the animal is moving inside the ring it would always be 'moving away from' the zone as the zone surrounds it, but it could still either be getting closer or further from the zone. (In fact, it would always be moving both towards and away from the zone, so because of this ambiguity, ANY-maze simply wouldn't score either of these measures.)

3.60 Initial heading error to the zone

<i>Description</i>	Reports the angle between the animal's heading at the start of the test and a direct heading to the zone.
<i>Calculation method</i>	<p>Calculation of this measure depends on the settings in the Analysis options > Heading error sub-element of the protocol. Specifically, the options control both how the animal's heading at the start of the test is determined and what part of the zone is used to calculate the heading error.</p> <p>There are two options for how the animal's heading at the start of the test is determined - one uses a specific time delay, the other a specific distance. In the first case, the animal's heading is taken to be the vector from its first position in the test to the first position detected after the specified time interval has elapsed. In the second case, the heading is taken to be the vector from the animal's first position in the test to the first position that's more than the specified distance from it. In both cases, positions that are detected while the animal is considered to be immobile (if immobility detection is switched on) are ignored - thus in the first case, the animal must be <i>mobile</i> for the period that is specified.</p> <p>Having determined the animal's initial heading, the system then calculates the heading from the first position in the test to the zone. To do this, ANY-maze can use one of two methods (again, these are specified using the Analysis options > Heading error sub-element in the protocol); it can either simply calculate the heading to the centre of the zone or it can calculate the heading to any part of the zone.</p> <p>In the first case (see figure 5), the <i>centre</i> of the zone is taken to be the zone's 'centre of mass' (i.e. the mean x, y coordinate of all the points in the zone) and the heading error is defined as the angle between this heading and the animal's initial heading.</p>




Figure 5. Calculation of the initial heading error using the centre of the zone: The animal's initial heading error is the angle 'A' between its initial heading and the direct heading to the centre of the zone.

It's important to understand that a zone's centre of mass may actually be outside the zone. For example, consider a ring-shaped zone; the centre of mass will be in the centre of the ring, but this point will not be within the zone.

The second method of calculating the heading to the zone is to consider the heading to every position on the zone's perimeter - in this case, the heading error is the smallest angle between the animal's heading and the heading to any perimeter point - see figure 6.




Figure 6. Calculation of the initial heading error using the entire area of the zone: The animal's initial heading error is the smallest angle 'A' between its initial heading and any point on the zone perimeter.

In the case of small zones the difference between the two calculation techniques is usually minimal, but for a large zone, it can make a substantial difference to the result.

Analysis across time This measure cannot be analysed across time.

<i>Units</i>	Degrees
<i>Notes</i>	None

3.61 Signed initial heading error to the zone

Description Reports the angle between the animal's heading at the start of the test and a direct heading to the zone. A positive initial heading error means the zone is to the animal's right, and a negative initial heading error signifies the zone is to the animal's left.

Calculation method Calculation of this measure depends on the settings in the Analysis options > Heading error sub-element of the protocol. Specifically, the options control both how the animal's heading at the start of the test is determined and what part of the zone is used to calculate the heading error.

There are two options for how the animal's heading at the start of the test is determined - one uses a specific time delay, the other a specific distance. In the first case, the animal's heading is taken to be the vector from its first position in the test to the first position detected after the specified time interval has elapsed. In the second case, the heading is taken to be the vector from the animal's first position in the test to the first position that's more than the specified distance from it. In both cases, positions that are detected while the animal is considered to be immobile (if immobility detection is switched on) are ignored - thus in the first case, the animal must be *mobile* for the period that is specified.

Having determined the animal's signed initial heading, the system then calculates the heading from the first position in the test to the zone. To do this, ANY-maze can use one of two methods (again, these are specified using the Analysis options > Heading error sub-element in the protocol); it can either simply calculate the heading to the centre of the zone, or it can calculate the heading to any part of the zone.

In the first case (see figure 7), the *centre* of the zone is taken to be the zone's 'centre of mass' (i.e. the mean *x, y* coordinate of all the points in the zone) and the heading error is defined as the angle between this heading and the animal's initial heading.




Figure 7. Calculation of the signed initial heading error using the centre of the zone: The animal's signed initial heading error is the angle 'A' between its initial heading and the direct heading to the centre of the zone. In this figure, the signed initial heading error will be a negative value, as the zone is to the animal's left.

It's important to understand that a zone's centre of mass may actually be outside the zone. For example, consider a ring-shaped zone; the centre of mass will be in the centre of the ring, but this point will not be within the zone.

The second method of calculating the heading to the zone is to consider the heading to every position on the zone's perimeter - in this case, the heading error is the smallest angle between the animal's heading and the heading to any perimeter point - see figure 8.




Figure 8. Calculation of the signed initial heading error using the entire area of the zone: The animal's initial heading error is the smallest angle 'A' between its initial heading and any point on the zone perimeter.

In the case of small zones, the difference between the two calculation techniques is usually minimal, but for a large zone it can make a substantial difference to the result.

<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	None

3.62 Average absolute heading error to the zone

<i>Description</i>	Reports the average absolute angle between the animal's heading and a direct heading to the zone.
<i>Calculation method</i>	The method used to calculate this measure depends on the option specified in the Analysis options > Heading error sub-element in the protocol. Specifically, there

are two ways to determine the heading to the zone - using the centre of the zone or using the entire zone area.

In the first case, the heading to the zone is taken to be the heading to the *centre* of the zone, where the centre is defined as the zone's 'centre of mass' (i.e. the mean x, y coordinate of all the points in the zone). It's important to understand that the centre of mass may actually be outside the zone. For example, consider a ring-shaped zone; the centre of mass will be in the centre of the ring, but this point will not be within the zone. With this definition of the heading to the zone, the heading error for a position is calculated as follows: The animal's heading is defined as the vector that joins the position with the next position in time. The heading to the zone is defined as the vector that joins the position to the centre of the zone and the heading error is defined as the angle between the two vectors - see figure 9.




Figure 9. Calculation of the average heading error using the centre of the zone: The animal's heading error is the angle 'A' between its heading and the direct heading to the centre of the zone.

In the case where the heading to the zone is defined using the entire zone area, the calculation of the heading error is performed as follows: The animal's heading is defined as the vector that joins the position with the next position in time. The heading to the zone is then calculated for every point on the zone's perimeter, and the angle between this heading and the animal's heading is calculated. The smallest angle is the heading error - see figure 10.




Figure 10. Calculation of the average heading error using the entire area of the zone: The animal's heading error is the smallest angle 'A' between its heading and the direct heading to any point on the zone perimeter.

Irrespective of which method is used to calculate the individual heading errors, the *average* absolute heading error is calculated in the same way: Each absolute heading error angle is multiplied by the time for which it persisted (i.e. the time from one position of the animal to the next). This product is then summed for the entire test (or time period). The final sum is then divided by the test duration (or the duration of the time period). This seemingly strange method of calculating the average is required because positions in ANY-maze are not necessarily recorded at a fixed frequency.

If immobility is being detected in a test, then all positions when the animal is deemed to be immobile are ignored in the calculation of the average heading error. If immobility is not being detected, then all positions are used with the caveat that a position must be at least a minimum distance from the previous position for it to be considered. The value used for this minimum distance is based on the size of the animal.

Analysis across time This measure can be analysed across time. The result is based on just those positions that fall within the time period.

Units Degrees.

Notes None

3.63 Time moving towards the zone

Description Reports the total amount of time for which the animal was moving towards the zone.

Calculation method The method used to calculate this measure depends on the option specified in the Orientation and movement element of the protocol. Specifically, the measure can either be based on the centre of the zone or on the zone's entire area.

In the first case, the measure is calculated as follows: For each position of the animal, a vector is created between the current position and the next position. A

second vector is then created between the current position and the centre of the zone. The angle between these two vectors is calculated - see figure 11.




Figure 11. ANY-maze calculates the angle 'A' between the animal's heading and a direct heading to the zone; if this is less than a critical angle, the animal is deemed to be moving towards the zone.

Here, the *centre* of the zone is defined as the zone's 'centre of mass' (i.e. the mean x, y coordinate of all the points in the zone). It's important to understand that the centre of mass may actually be outside the zone. For example, consider a ring-shaped zone; the centre of mass will be in the centre of the ring, but this point will not be within the zone.

Having calculated the angle between the two vectors ('A'), it is compared to the critical angle for movement towards a zone (see notes). If the angle is less than this critical angle, then the animal is deemed to be moving towards the zone (unless it is *also* moving away from it - see notes) and the time it took to move from the current position to the next position is added to the total time moving towards the zone.

The second method used to calculate this measure uses the entire area of the zone. In this case, for each position of the animal, a vector is created between the current position and the next position. All possible vectors from the current position to the points on the zone's perimeter are then calculated, and the angle between each one and the animal's heading vector is calculated. The smallest of these angles is found. This angle is then compared to the critical angle for movement towards a zone in the same way as for the first calculation method (see above).

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving towards the zone is specified in the protocol's Orientation and movement element. In fact, the value entered in the analysis options is twice the critical angle described here, as this is more intuitive. The default critical angle is 90° (i.e. by default, the angle 'A' will be compared to 45°).

When calculating whether the animal is moving towards a zone, the system takes into consideration whether the animal is *also* moving away from the zone (see Time moving away from the zone). If it is, then the position is deemed to be ambiguous and the time from the position to the next position is NOT added to the overall result. A simple example of when this situation would arise is with a ring-shaped zone. If the animal is inside the ring, then no matter what direction it moves in, it will be moving both towards the zone and away from it (as the zone surrounds the animal). The measures Time getting closer to the zone and Time getting further away from the zone offer alternatives that avoid this problem.

3.64 Time moving away from the zone

Description Reports the total amount of time for which the animal was moving away from the zone.

Calculation method The method used to calculate this measure depends on the option specified in the Orientation and movement element of the protocol. Specifically, the measure can either be based on the centre of the zone or on the zone's entire area.

In the first case, the measure is calculated as follows: For each position of the animal, a vector is created between the current position and the next position. A second vector is then created between the current position and the centre of the zone. The angle between these two vectors ('B') is calculated and deducted from 180°, yielding angle 'A' - see figure 12.




Figure 12. ANY-maze calculates the angle 'B' between the animal's heading and a direct heading to the zone; this is deducted from 180°, yielding angle 'A'. If this angle is less than a critical angle, the animal is deemed to be moving away from the zone.

Here, the *centre* of the zone is defined as the zone's 'centre of mass' (i.e. the mean x, y coordinate of all the points in the zone). It's important to understand that the centre of mass may actually be outside the zone. For example, consider a ring-

shaped zone; the centre of mass will be in the centre of the ring, but this point will not be within the zone.

Having calculated the angle 'A', it is compared to the critical angle for movement away from a zone (see notes). If the angle is less than this critical angle, then the animal is deemed to be moving away from the zone (unless it is *also* moving towards it - see notes) and the time it took to move from the current position to the next position is added to the total time moving away from the zone.

The second method used to calculate this measure uses the entire area of the zone. In this case, for each position of the animal, a vector is created between the current position and the next position. All possible vectors from the current position to the points on the zone's perimeter are then calculated, and the angle between each one and the animal's heading vector is calculated. The largest of these angles is found. This angle is deducted from 180° and the result is compared to the critical angle for movement away from a zone in the same way as for the first calculation method (see above).

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving away from the zone is also specified in the protocol's Orientation and movement element. In fact, the value entered in the analysis options is twice the critical angle as described here, as this is more intuitive. The default critical angle is 90° (i.e. by default, the angle 'A' will be compared to 45°).

When calculating whether the animal is moving away from a zone, the system takes into consideration whether the animal is *also* moving towards the zone (see Time moving towards the zone). If it is, then the position is deemed to be ambiguous and the time from the position to the next position is NOT added to the overall result. A simple example of when this situation would arise is a ring-shaped zone. If the animal is inside the ring, then no matter what direction it moves it will be moving both towards the zone and away from it (as the zone surrounds the animal). The measures Time getting closer to the zone and Time getting further away from the zone offer alternatives that avoid this problem.

3.65 Time oriented towards the zone

Description Reports the total amount of time the animal was oriented towards the zone (while outside the zone).

Calculation method The animal's orientation is defined by a vector connecting the *start of the head vector position* to the head position. For most animals, excluding quails, the *start of the head vector position* is the animal's centre. However, in the case of quails, it is a point near the base of the neck, which was determined during tracking. A second vector from the animal's head position to a point on the border of the zone is determined, and if the absolute angle between the vectors is less than or equal to a *critical angle* (see notes), then the animal is deemed to be oriented towards the zone. This is repeated for every point on the zone's border until either the animal is

found to be oriented towards the zone or all the points on the zone's border have been processed.

If the animal is oriented towards the zone then the time between the previous position of the animal and its current position is added to the total time oriented towards the zone (irrespective of whether it was oriented towards the zone at the previous position).

Analysis across time This measure can be analysed across time.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving towards the zone is also specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle as described here, as this is more intuitive. The default critical angle is 60° (i.e. by default, the angle between the vectors will be compared to 30°).

3.66 Time oriented towards the centre of the zone when inside zone

Description Reports the amount of time the animal was oriented towards the centre of the zone while it was inside the zone.

Calculation method The animal's orientation is defined by a vector connecting the *start of the head vector position* to the head position. For most animals, excluding quails, the *start of the head vector position* is the animal's centre. However, in the case of quails, it is a point near the base of the neck, which was determined during tracking. A second vector from the animal's head position to the centre of the zone is determined, and if the angle between the vectors is less than or equal to a *critical angle* (see notes), then the animal is deemed to be oriented towards the zone. The amount of time for which this is the case while the animal is inside the zone (as determined by the zone entry criteria) is summed.

Analysis across time This measure can be analysed across time.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving towards the zone is also specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle as described here, as this is more intuitive. The default critical angle is 60° (i.e. by default, the angle between the vectors will be compared to 30°).

3.67 Absolute turn angle while in the zone

Description Reports the sum of the absolute angle between each movement vector of the animal while it was inside the zone.

Calculation method For each position of the animal that is inside the zone, a vector of movement from one position of the animal's centre point to the next is created. For each vector, the angle between it and the previous vector is calculated with anti-clockwise movement being negative and clockwise movement being positive (i.e. the angle is

from -180° to 180°). The absolute value of this angle is summed for all the positions of the animal within the zone throughout the test or time period.

<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions within the specific time period.
<i>Units</i>	Degrees
<i>Notes</i>	From this measure, it is easy to use calculations to derive measures such as <i>Meander in the zone</i> and <i>Angular velocity in the zone</i> . The former is the <i>Absolute turn angle while in the zone</i> divided by the <i>Distance travelled in the zone</i> and the latter is the <i>Absolute turn angle while in the zone</i> divided by the <i>Time in the zone</i> .

3.68 Absolute head turn angle while in the zone

<i>Description</i>	Reports the cumulative absolute angle through which the animal's head moved while it was in the zone. For example, if while in the zone, the animal moved its head 30° to the left and then moved its head 45° to the right, the absolute head turn angle would be 75°.
<i>Calculation method</i>	For each position of the animal's head, a vector is created from the animal's centre point to the head. The angle between this vector and the same vector for the previous position of the animal's head is calculated, and the absolute value of this angle is summed whenever the animal is in the zone.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	<p>As described above, calculation of this measure requires two vectors which are generated from two consecutive head positions. It is the second of these two positions which is used to determine whether the animal is in the zone. Thus, for example, the animal might be oriented North and standing outside the zone; if it then moves into the zone and is then oriented West, the 90° change in orientation will all be attributed to the zone that it has just entered.</p> <p>Whether the animal is in the zone is determined by the zone entry settings. Note that these may not require that either the head or the centre point are actually in the zone.</p>

3.69 Freezing bouts in the zone

<i>Description</i>	Reports the number of times the animal <i>froze</i> while in the zone.
<i>Calculation method</i>	Each time the animal starts to freeze, a check is made to determine whether it is in the zone. If it is, the count is incremented by one.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	None

3.70 Time freezing in the zone

<i>Description</i>	Reports the total time the animal was <i>frozen</i> while in the zone.
<i>Calculation method</i>	The duration of each bout of freezing is calculated. If the animal was in the zone while frozen, then the duration is added to the result for the zone.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

3.71 Freezing bouts while investigating the zone

<i>Description</i>	Reports the number of times the animal <i>froze</i> while investigating the zone.
<i>Calculation method</i>	Each time the animal starts to freeze, a check is made to determine whether it is investigating the zone. If it is, the count is incremented by one.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.72 Time freezing while investigating the zone

<i>Description</i>	Reports the total time that the animal was <i>frozen</i> while investigating the zone.
<i>Calculation method</i>	The duration of each bout of freezing is calculated. If the animal was investigating the zone while frozen, then the duration is added to the result for the zone.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the zone is an investigation zone.

3.73 Number of rears in the zone

<i>Description</i>	Reports the number of times the animal reared while in the zone.
<i>Calculation method</i>	Depends on the method used to detect zone entries (and thus by implication, zone exits too) - see Choosing how ANY-maze should detect entries into a zone for more details.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

3.74 Total time rearing in the zone

<i>Description</i>	Reports the total amount of time for which the animal was rearing while it was in the zone.
<i>Calculation method</i>	Sums the duration of each bout of rearing that occurred while the animal was in the zone. If the animal enters the zone when it is already rearing, then the time will be counted from the time the animal entered the zone and not when the rearing bout started. If the animal exits the zone while rearing, then the time will stop at the time the animal exits the zone and not at the end of the rearing bout. For these reasons, it is possible for the result of this measure to be non-zero when the result for Number of rears in the zone is zero.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

3.75 Latency to first rear in the zone

<i>Description</i>	Reports the latency to first time that the animal reared in the zone.
<i>Calculation method</i>	The time when the first bout of rearing started while the animal was in the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	<p>This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.</p> <p>This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis options element.</p>

3.76 Average duration of a rear in the zone

<i>Description</i>	Reports the average duration of the rearing bouts in the zone.
<i>Calculation method</i>	The result of Total time rearing in the zone divided by Number of rears in the zone.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	<p>This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.</p> <p>This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.</p>

3.77 Maximum duration of a rear in the zone

<i>Description</i>	Reports the duration of the longest bout of rearing in the zone.
<i>Calculation method</i>	The duration of each bout of rearing in the zone is calculated and the longest bout is found. Note that a bout of rearing in the zone starts when the animal is in the zone and begins to rear OR when the animal enters the zone when it is already rearing. Similarly, a bout ends when the animal is in the zone and stops rearing OR the animal exits the zone while rearing.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

3.78 Minimum duration of a rear in the zone

<i>Description</i>	Reports the duration of the shortest bout of rearing in the zone.
<i>Calculation method</i>	The duration of each bout of rearing in the zone is calculated and the shortest bout is found. Note that a bout of rearing in the zone starts when the animal is in the zone and begins to rear OR when the animal enters the zone when it is already rearing. Similarly, a bout ends when the animal is in the zone and stops rearing OR the animal exits the zone while rearing.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the apparatus is being viewed from the side. ANY-maze actually detects rearing by analysing the shape of the animal, and therefore this measure will only work reliably if there is good contrast between the animal and the background of the apparatus.

3.79 Time spent in Whishaw's Corridor

<i>Description</i>	Reports the amount of time the animal spent in the Whishaw's Corridor for the zone.
<i>Calculation method</i>	The Whishaw's Corridor for the zone is determined based on the start position of the animal in the test, the centre point of the zone, and the width of the corridor specified in the protocol - see figure 13.




Figure 13. Example of the Whishaw's Corridor in a water-maze. The corridor is centred on a line (shown in red) running from the animal's start position to the centre of the platform zone. The corridor itself (shown in blue) has a width specified in the protocol.

The time spent in the corridor is then calculated by summing the duration of each visit to the corridor, where a visit starts at the time of entry and ends at the time of exit.

<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the total amount of time the animal spent within the corridor during that time period.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is only available if the Whishaw's Corridor width has been specified for the zone - see setting up a zone.

3.80 Distance travelled in Whishaw's Corridor

<i>Description</i>	Reports the distance travelled by the animal in the Whishaw's Corridor for the zone.
<i>Calculation method</i>	The Whishaw's Corridor for the zone is determined based on the start position of the animal in the test, the centre point of the zone, and the width of the corridor specified in the protocol - see figure 14.




Figure 14. Example of Whishaw's Corridor in a water-maze. The corridor is centred on a line (shown in red) running from the animal's start position to the centre of the platform zone. The corridor itself (shown in blue) has a width specified in the protocol.

The distance travelled in the corridor is then calculated by summing the distance travelled during each visit to the corridor, where a visit starts when an animal enters the zone and ends when it exits the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the distance travelled within the corridor during that time period.

Units Metres

Notes This measure is only available if the Whishaw's Corridor width has been specified for the zone - see Setting up a zone.

In common with the way in which distance travelled in a zone is calculated, the distance travelled by the animal between a point outside the corridor and a point inside the corridor (i.e. when entering the corridor) is NOT included in the distance travelled in the corridor, whereas the distance travelled by the animal between a point inside the corridor and a point outside the corridor (i.e. when exiting the corridor) IS included in the distance travelled in the corridor.

3.81 Time the animal's head was in the zone when its centre was outside the zone

Description Reports the amount of time for which the animal's head position was in the zone while its centre position was outside the zone

Calculation method For every position of the animal's head, a check is made to determine whether the head position is in the zone and the centre position is outside the zone (note that this determination does not use the zone entry criteria). The duration of each occurrence is calculated, and the total duration of all occurrences is summed.

<i>Analysis across time</i>	This measure can be analysed across time
<i>Units</i>	Seconds
<i>Notes</i>	None

3.82 Path efficiency to first entry to the zone

<i>Description</i>	This measure represents an index of the efficiency of the path taken by the animal to get from the first position in the test to the first position within the zone. A value of 1 indicates perfect efficiency - the animal moved in a straight line - values less than 1 indicate decreasing efficiency.
<i>Calculation method</i>	The straight line distance between the first position in the test and the first position in the zone is divided by the total distance travelled by the animal until it first entered the zone
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	None
<i>Notes</i>	<p>This measure is intended for use in water-maze experiments but is available in all tests.</p> <p>This result is affected by the setting in Analysis option to Use zero as the result for undefined path efficiency to zones.</p> <p>This measure will be undefined, and therefore reported as #N/A (or zero, see previous note), if the animal's route to the first entry to this zone passed through a hidden zone. This is because ANY-maze can't know the distance the animal travelled while in the hidden zone and therefore can't calculate the path efficiency.</p>

3.83 Corrected integrated path length

<i>Description</i>	Reflects how efficiently the animal moved from the start position to the zone, lower scores are better.
<i>Reference</i>	Barnes CA, et al. (1997) Multistability of cognitive maps in the hippocampus of old rats <i>Nature</i> 388 : 272-5
<i>Calculation method</i>	The difference between the sum of the sampled distances from the target zone and the shortest possible sum if the animal had moved directly to the zone at its mean speed (where the mean speed is the mean speed prior to entering the zone). The reported value applies to the first entry to the zone (see notes).
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Metres · seconds
<i>Notes</i>	By default, this value is based on the animal's path until the first zone entry, which means that if the animal never enters the zone the result will be reported as #N/A. This can be altered by selecting the Analysis option to Calculate CIPL for a zone based on the animal's path throughout the test, which calculates a zone's CIPL irrespective of how many times the animal entered the zone - follow the link for more details.

This value can be negative. Consider a situation in which the animal starts 10m from the target zone and moves to it in a direct line in a period of 10s - the integrated path length (i.e. the area under the curve of distance vs. time) would be 50m.s. But if the animal had actually moved 9m in the first 1s, so it was then 1m from the zone, and then moved the remaining 1m over the following 9s, then the actual integrated path length would be 10m.s and therefore the CIPL would yield a result of $10 - 50 = -40\text{m.s}$

3.84 Number of line crossings while in the zone

<i>Description</i>	Reports the number of times the animal's centre point moved from one area of the apparatus map to another - i.e. crossed the lines which constitute the map, while the animal was in the specific zone.
<i>Calculation method</i>	The apparatus is divided into unique areas by the apparatus map. For each animal position recorded in the experiment, the area which contains the animal's centre point is found. Each time this changes, if the animal is in the zone, the measure's value is increased by 1.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	<p>When an animal enters a zone, it will usually cross a line; these line crossings are counted in the zone the animal has <i>entered</i>. This means that if the entire apparatus is divided into zones and the number of line crossings for all the zones is summed, it will equal the total number of line crossings for the apparatus as a whole. If you want to know the number of line crossings <i>within</i> the zone (i.e. excluding the line crossings that occur when the animal enters the zone), then use a calculation to subtract the number of zone entries from the number of line crossings in the zone.</p> <p>It's important to understand that ANY-maze uses the animal's centre point when calculating the measure and therefore it can be prone to 'spurious entries' if an animal straddles a line between two areas (i.e. by moving a very small amount, the animal can apparently cross a line many times). This problem can be overcome by setting one zone for each area, using the percentage of the animal that's in the zone to score zone entries, and then using a calculation to sum all the entries into these zones.</p>

4 Point measures

4.1 Average distance from the point

<i>Description</i>	Reports the average distance from the animal to the point.
<i>Calculation method</i>	For each position of the animal, the distance is calculated from the point to the animal; this distance is averaged throughout the entire test or time period. The method actually used to calculate this distance depends on the setting made when setting up the point. Specifically, the distance can either be calculated based on the part of the animal which is closest to the point (excluding its tail), or based on the position of the centre of the animal.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the average distance from the point to the animal during the time period.
<i>Units</i>	Metres
<i>Notes</i>	None

4.2 Maximum distance from the point

<i>Description</i>	Reports the maximum distance from the animal to the point.
<i>Calculation method</i>	For each position of the animal, the distance is calculated from the point to the part of the animal that's closest to the point (the animal's tail is excluded). The maximum distance during the entire test or time period is the result. The method actually used to calculate this distance depends on the setting made when setting up the point. Specifically, the distance can either be calculated based on the part of the animal which is closest to the point (excluding its tail), or based on the position of the centre of the animal.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the maximum distance from the animal to the point during the time period.
<i>Units</i>	Metres
<i>Notes</i>	None

4.3 Minimum distance from the point

<i>Description</i>	Reports the minimum distance from the animal to the point.
<i>Calculation method</i>	For each position of the animal, the distance is calculated from the point to the part of the animal that's closest to the point (the animal's tail is excluded). The minimum distance during the entire test or time period is the result. The method actually used to calculate this distance depends on the setting made when setting up the point. Specifically, the distance can either be calculated based on the part of the animal which is closest to the point (excluding its tail), or based on the position of the centre of the animal.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the minimum distance from the animal to the point during the time period.

<i>Units</i>	Metres
<i>Notes</i>	None

4.4 Average distance of the animal's head from the point

<i>Description</i>	Reports the average distance from the animal's head to the point.
<i>Calculation method</i>	Calculates the distance from the point to the position of the animal's head. Averages this distance through the entire test or time period.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the average distance from the point to the animal's head during the time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if Head tracking is turned on.

4.5 Maximum distance of the animal's head from the point

<i>Description</i>	Reports the maximum distance from the animal's head to the point.
<i>Calculation method</i>	Calculates the distance from the point to the position of the animal's head. The maximum distance during the entire test or time period is the result.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the maximum distance from the animal's head to the point during the time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if Head tracking is turned on.

4.6 Minimum distance of the animal's head from the point

<i>Description</i>	Reports the minimum distance from the animal's head to the point.
<i>Calculation method</i>	Calculates the distance from the point to the position of the animal's head. The minimum distance during the entire test or time period is the result.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the minimum distance from the animal's head to the point during the time period.
<i>Units</i>	Metres
<i>Notes</i>	This measure is only available if Head tracking is turned on.

4.7 Time moving towards the point

<i>Description</i>	Reports the total amount of time for which the animal was moving towards the point.
<i>Calculation method</i>	For each position of the animal, a vector is created between the current position and the next position. A second vector is then created between the current position and the point. The angle between these two vectors is calculated - see figure 15.




Figure 15. ANY-maze calculates the angle 'A' between the animal's heading and a direct heading to the point.

The absolute angle between the two vectors, 'A', is compared to the critical angle for movement towards a point (see notes). If the angle is less than this critical angle, then the animal is deemed to be moving towards the point. If the animal is moving towards the point, then the time taken to move from the current position to the next position is added to the total time moving towards the point.

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving towards the point is specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 90° (i.e. by default the angle 'A' will be compared to 45°).

4.8 Time moving away from the point

Description Reports the total amount of time for which the animal was moving away from the point.

Calculation method For each position of the animal, a vector is created between the current position and the next position. A second vector is then created between the current position and the point. The angle between these two vectors is calculated and subtracted from 180° - see figure 16.




Figure 16. ANY-maze calculates the angle 'B' between the animal's heading and a direct heading to the point. This is then subtracted from 180°, yielding angle 'A'.

The angle 'A' is compared to the critical angle for movement away from a point (see notes). If the value is less than this critical angle, then the animal is deemed to be moving away from the point. If the animal is moving away from the point, then the time taken to move from the current position to the next position is added to the total time moving away from the point.

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving away from the point is specified in the protocol's Orientation and movement element. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 90° (i.e. by default the angle 'A' will be compared to 45°).

4.9 Time the animal's head was moving towards the point

Description Reports the total amount of time for which the animal's head was moving towards the point.

Calculation method For each position of the animal's head, a vector is created between the current position and the next position. A second vector is then created between the current position and the point. The angle between these two vectors is calculated - see figure 17.




Figure 17. ANY-maze calculates the angle 'A' between the heading of the animal's head and the direct heading to the point.

The angle 'A' is compared to the critical angle for movement towards a point (see notes). If the angle is less than this critical angle, then the animal's head is deemed to be moving towards the point. If the animal's head is moving towards the point, then the time it took to move from the previous position to the current position is added to the total time moving towards the point.

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is moving towards the point is specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 90° (i.e. by default the angle *beta* is compared to 45°).

This measure is only available if Head tracking is turned on.

4.10 Time the animal's head was moving away from the point

Description Reports the total amount of time for which the animal's head was moving away from the point.

Calculation method For each position of the animal's head, a vector is created between the current position and the next position. A second vector is then created between the current position and the point. The angle between these two vectors is calculated and the result is subtracted from 180° - see figure 18.




Figure 18. ANY-maze calculates the angle 'B' between the heading of the animal's head and the direct heading to the point; this is then subtracted from 180°, yielding angle 'A'.

The angle 'A' is compared to the critical angle for movement away from a point (see notes). If the angle is less than this critical angle, then the animal's head is deemed to be moving away from the point. If the animal's head is moving away from the point, then the time it took to move from the previous position to the current position is added to the total time moving away from the point.

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal's head is moving away from the point is specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 90° (i.e. by default the angle 'A' will be compared to 45°).

This measure is only available if Head tracking is turned on.

4.11 Average speed moving towards the point

Description Reports the average speed of the animal when it was moving towards the point.

Calculation method For each position of the animal, a vector is created between the current position and the next position. A second vector is then created between the current position and the point. The angle between these two vectors is calculated - see figure 19.




Figure 19. ANY-maze calculates the angle 'A' between the animal's heading and a direct heading to the point.

The absolute angle between the two vectors, 'A', is compared to the critical angle for movement towards a point (see notes). If the angle is less than this critical angle, then the animal is deemed to be moving towards the point. If the animal is moving towards the point, then the time taken to move from the current position to the next position is added to the total time moving towards the point and the distance from the current position to the next position is added to the total distance moving towards the point. At the end of the test the total distance moving towards the point is divided by the total time moving towards the point, yielding the average speed moving towards the point.

<i>Analysis across time</i>	This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.
<i>Units</i>	Metres per second
<i>Notes</i>	The <i>critical angle</i> used to define whether the animal is moving towards the point is specified in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 90° (i.e. by default the angle 'A' will be compared to 45°).

4.12 Initial heading error to the point

<i>Description</i>	Reports the angle between the animal's heading at the start of the test and a direct heading to the point.
<i>Calculation method</i>	Calculation of this measure depends on the settings in the Analysis options > Heading error sub-element of the protocol. Specifically, there are two options for how the animal's heading at the start of the test is determined - one uses a specific time delay, the other a specific distance. In the first case, the animal's heading is taken to be the vector from its first position in the test to the first position detected after the specified time interval has elapsed. In the second case, the heading is taken to be the vector from the animal's first position in the test to the first position that's more than the specified distance from it. In both cases, positions that are detected while the animal is considered to be immobile (if immobility

detection is switched on) are ignored - thus in the first case, the animal must be *mobile* for the period that is specified.

Having determined the animal's initial heading, the system then calculates the vector from the first position in the test to the point. The angle between this vector and the animal's heading vector is the initial heading error.

<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Degrees
<i>Notes</i>	None

4.13 Average absolute heading error to the point

<i>Description</i>	Reports the average absolute angle between the animal's heading and a direct heading to the point.
<i>Calculation method</i>	<p>The animal's heading is defined as the vector that joins the position with the next position in time. The heading to the point is defined as the vector that joins the animal's position to the point, and the heading error is defined as the angle between the two vectors. This angle is calculated for every position of the animal, and the angle's absolute value is summed and then divided by the number of positions.</p> <p>If immobility is being detected in a test, then all positions when the animal is deemed to be immobile are ignored in the calculation of the average heading error. If immobility is not being detected, then all positions are used - with the caveat that a position must be at least a minimum distance from the previous position for it to be considered. The minimum distance used is based on the size of the animal.</p>
<i>Analysis across time</i>	This measure can be analysed across time. The result is based on just those positions that fall within the time period.
<i>Units</i>	Degrees
<i>Notes</i>	None

4.14 Time the animal's head was oriented towards the point

<i>Description</i>	Reports the total amount of time for which the animal's head was oriented towards the point.
<i>Calculation method</i>	The animal's orientation is defined by a vector connecting the <i>start of the head vector position</i> to the head position. For most animals, excluding quails, the <i>start of the head vector position</i> is the animal's centre. However, in the case of quails, it is a point near the base of the neck, which was determined during tracking. A second vector is created between the position of the animal's head and the position of the point. The angle between these two vectors is calculated - see figure 20.




Figure 20. ANY-maze calculates the angle beta between the two vectors V_a (the vector which defines the animal's orientation) and V_b (the vector from the animal to the point).

The absolute angle between the two vectors (*beta*) is compared to the critical angle for movement away from/towards a point (see notes). If the angle is less than this critical angle, then the animal is deemed to be oriented towards the point. If the animal is oriented towards the point, then the time between the previous position of the animal and its current position is added to the total time oriented towards the point (irrespective of whether it was oriented towards the point at the previous position).

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal is oriented towards the point is the same angle as is used to determine whether the animal is moving towards or away from a point, and is defined in the Orientation and movement element of the protocol. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 60° (i.e. by default the angle *beta* will be compared to 30°).

This measure is only available if Head tracking is turned on.

4.15 Time the animal's head was oriented away from the point

Description Reports the total amount of time for which the animal's head was oriented away from the point.

Calculation method The animal's orientation is defined by a vector connecting the *start of the head vector position* to the head position. For most animals, excluding quails, the *start of the head vector position* is the animal's centre. However, in the case of quails, it is a point near the base of the neck, which was determined during tracking. A second vector is created between the position of the animal's head and the position of the point. The angle between these two vectors is calculated - see figure 21.




Figure 21. ANY-maze calculates the angle beta between the two vectors V_a (the vector which defines the animal's orientation) and V_b (the vector from the animal to the point).

The angle between the two vectors (*beta*) is subtracted from 180° and the absolute value of the result is calculated; this is compared to the critical angle for movement towards/away from a point (see notes). If the angle is less than this critical angle, then the animal is deemed to be oriented away from the point. If the animal is oriented away from the point, then the time between the previous position of the animal and its current position is added to the total time oriented away from the point (irrespective of whether it was orientated away from the point at the previous position).

Analysis across time This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.

Units Seconds

Notes The *critical angle* used to define whether the animal's head is oriented away from the point is the same angle as is used to determine whether the animal is moving towards or away from a point, and is defined in the protocol's Orientation and movement element. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 60° (i.e. by default the angle $abs(180-beta)$ will be compared to 30°).

This measure is only available if Head tracking is turned on.

4.16 Number of times the animal's head was oriented towards the point

Description Reports the count of occasions when the animal's head was oriented towards the point.

Calculation method The animal's orientation is defined by a vector connecting the *start of the head vector position* to the head position. For most animals, excluding quails, the *start of the head vector position* is the animal's centre. However, in the case of quails, it is a point near the base of the neck, which was determined during tracking. A second

vector is created between the position of the animal's head and the position of the point. The angle between these two vectors is calculated - see figure 22.




Figure 22. ANY-maze calculates the angle beta between the two vectors V_a (the vector which defines the animal's orientation) and V_b (the vector from the animal to the point).

The absolute angle between the two vectors (*beta*) is compared to the critical angle for movement away from/towards a point (see notes). If the angle is less than this critical angle, then the animal is deemed to be oriented towards the point. If the animal is now oriented towards the point when previously it wasn't, then the count of occasions when the animal's head was oriented towards the point is increased by one.

<i>Analysis across time</i>	This measure can be analysed across time. The result is calculated using just those positions which fall within the time period.
<i>Units</i>	Seconds
<i>Notes</i>	The <i>critical angle</i> used to define whether the animal is oriented towards the point is the same angle as is used to determine whether the animal is moving towards or away from a point, and is defined in the protocol's Orientation and movement element. In fact, the value entered in the analysis options is twice the critical angle, as this is more intuitive. The default critical angle is 60° (i.e. by default the angle <i>beta</i> will be compared to 30°). This measure is only available if Head tracking is turned on.

4.17 Point X coordinate

<i>Description</i>	Reports the X coordinate of the location of a point whose position is determined by where the animal spent the longest time either located or doing something, such as freezing.
<i>Calculation method</i>	A heat map for whatever the point is based on is created and the position of the 'hottest' value in the heat map is determined.
<i>Analysis across time</i>	This measure cannot be analysed across time.

<i>Units</i>	Pixels
<i>Notes</i>	This measure is not available for points whose location in the apparatus is set within the protocol. This is because the location will be the same for all tests, and so reporting it as a measure is meaningless. If you want to know the location of such a point then it is reported in the Protocol report.

4.18 Point Y coordinate

<i>Description</i>	Reports the Y coordinate of the location of a point whose position is determined by where the animal spent the longest time either located or doing something such as freezing.
<i>Calculation method</i>	A heat map for whatever the point is based on is created and the position of the 'hottest' value in the heat map is determined.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Pixels
<i>Notes</i>	This measure is not available for points whose location in the apparatus is set within the protocol. This is because the location will be the same for all tests, and so reporting it as a measure is meaningless. If you want to know the location of such a point then it is reported in the Protocol report.

4.19 Approximate time at the point

<i>Description</i>	Reports the approximate time the animal spent at the location of a point whose position is determined by where the animal spent the longest time either located or doing something, such as freezing.
<i>Calculation method</i>	A heat map for whatever the point is based on is created and the time the animal spent at the 'hottest' value in the heat map is determined.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	The value reported is approximate because the way a heat map is generated is itself approximate - heat maps being intended to provide a visual, rather than analytical, representation of the animal's behaviour during a test. This measure is not available for points whose location in the apparatus is set within the protocol.

5 Sequence measures

5.1 Number of sequences

<i>Description</i>	Reports the number of times the sequence was performed.
<i>Calculation method</i>	Each time a sequence ends, the number of sequences is incremented by one.
<i>Analysis across time</i>	This measure can be analysed across time. A sequence is considered to occur in the time period in which it ends . For example, if a sequence starts at time 15 seconds and continues until time 45 seconds, then (for 30 second time periods) it would be counted in the time period 30-60 seconds because this is the period in which it ended.
<i>Units</i>	None
<i>Notes</i>	None

5.2 Total time performing sequences

<i>Description</i>	Reports the total amount of time that the animal was performing sequences during the test.
<i>Calculation method</i>	Sums the time taken to complete each sequence.
<i>Analysis across time</i>	This measure can be analysed across time. A sequence is considered to occur in the time period in which it ends , therefore the time performing the sequence is ALL attributed to the time period in which it ends. In an extreme case, this can yield a result in which the time performing a sequence in a time period is <i>longer</i> than the period itself. For example, if a sequence started at time 5 seconds and ended at time 55 seconds, then (for 30 second time periods), its duration - 50 seconds - would be counted in the time period 30-60 seconds because this is the period in which it ended, but the duration would actually be longer than the time period. This apparent strange way of calculating results for periods is done to ensure that the result for Average time to complete the sequence is correct.
<i>Units</i>	Seconds
<i>Notes</i>	None

5.3 Latency to start of first sequence

<i>Description</i>	Reports the length of time which elapsed before the first sequence in the test started.
<i>Calculation method</i>	Time at which first sequence starts.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis element.

5.4 Latency to completion of first sequence

<i>Description</i>	Reports the length of time which elapsed before the first sequence in the test ended.
<i>Calculation method</i>	Time at which the sequence first ended.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the option to <i>Use the test duration as the latency for events which don't occur</i> in the Analysis element.

5.5 Average time to complete the sequence

<i>Description</i>	Reports the average amount of time taken to complete the sequences.
<i>Calculation method</i>	Calculated by dividing the Total time performing sequences by the Number of sequences. If there are no sequences in the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result is Total time performing sequences in the period divided by the Number of sequences in the period. If no sequences were performed during a period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

5.6 Maximum duration of a sequence

<i>Description</i>	Reports the maximum amount of time taken to complete a single sequence.
<i>Calculation method</i>	The duration of each sequence is calculated when the sequence ends and the longest one is found. If no sequences were performed during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the duration of the longest sequence which ended during the time period. If no sequence ended during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	None

5.7 Minimum duration of a sequence

<i>Description</i>	Reports the minimum amount of time taken to complete a single sequence.
<i>Calculation method</i>	The duration of each sequence is calculated when the sequence ends and the shortest one is found. If no sequences were performed during the test, then the result is undefined.

<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the duration of the shortest sequence which ended during the time period. If no sequence ended during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	None

5.8 Total distance travelled during sequences

<i>Description</i>	Reports the total distance travelled by the animal while performing sequences in the test.
<i>Calculation method</i>	Sums the distance travelled during each sequence in the test.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the sum of the distance travelled for all the sequences that ended during the time period.
<i>Units</i>	Metres
<i>Notes</i>	None

5.9 Average distance travelled per sequence

<i>Description</i>	Reports the average distance travelled by the animal while performing each of the sequences in the test.
<i>Calculation method</i>	Calculated by dividing the Total distance travelled during sequences by the Number of sequences. If there were no sequences performed in the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the Total distance travelled during sequences which ended during the time period divided by the Number of sequences which ended during the time period. If there were no sequences which ended during the time period, then the result is undefined.
<i>Units</i>	Metres
<i>Notes</i>	This measure is affected by the option to <i>Use zero as the result for undefined averages</i> in the Analysis options element.

5.10 Maximum distance travelled during a sequence

<i>Description</i>	Reports the maximum distance travelled by the animal while performing a single sequence.
<i>Calculation method</i>	The distance travelled while performing each sequence is calculated and the greatest distance is found. If no sequences were performed in the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the greatest distance travelled during any sequence which ended during the time period. If no sequence ended during the time period, then the result is undefined.
<i>Units</i>	Metres

Notes None

5.11 Minimum distance travelled during a sequence

Description Reports the minimum distance travelled by the animal while performing a single sequence.

Calculation method The distance travelled while performing each sequence is calculated and the smallest distance is found. If no sequences were performed in the test, then the result is undefined.

Analysis across time This measure can be analysed across time. The result for a time period is the smallest distance travelled during any sequence which **ended** during the time period. If no sequence ended during the time period, then the result is undefined.

Units Metres

Notes None

5.12 Average speed during the sequence

Description Reports the average speed of the animal while performing the sequences in the test.

Calculation method Calculated by dividing Total distance travelled during sequences by the Total time performing sequences. If no sequences were performed in the test, then the result is undefined.

Analysis across time This measure can be analysed across time. The result for a time period is the Total distance travelled during sequences which **ended** during the period divided by the Total time performing sequences which **ended** during the period. If no sequence ended during the time period, then the result is undefined.

Units Metres per second

Notes This measure is affected by the option to *Use zero as the result for undefined averages* in the Analysis options element.

6 Key measures

Each measure is available for the apparatus as a whole (i.e. irrespective of where the behaviour occurred) and also for each defined zone.

6.1 Number of presses

<i>Description</i>	Reports the number of times the key was pressed down, i.e. the number of times the animal started to exhibit the key's behaviour.
<i>Calculation method</i>	Counts the number of key presses.
<i>Analysis in zones</i>	Counts the number of times the key was pressed when the animal was in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of times the key was pressed down during the time period.
<i>Units</i>	None
<i>Notes</i>	None

6.2 Time pressed

<i>Description</i>	Reports the total amount of time the key was pressed down, i.e. the total amount of time that the animal was exhibiting the key's behaviour.
<i>Calculation method</i>	Sums the amount of time for which the key was pressed.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the key was pressed in the zone (or for an investigation zone, while the animal was investigating the zone).</p> <p>For a particular zone, it's possible for the <i>Time pressed</i> to be non-zero while the Number of presses is zero. This can occur if the animal enters the zone when the key is pressed. In this case, the time the key is pressed will be registered but the key press itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the key was pressed during the period.</p> <p>For a particular time period, it's possible for the <i>Time pressed</i> during the period to be non-zero while the Number of presses for the period is zero. This can occur if the key is already pressed at the start of the period. In this case, the time the key is pressed will be registered but the key press itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

6.3 Latency to first press

<i>Description</i>	Reports the amount of time that elapsed in the test before the key was pressed for the first time, i.e. the amount of time that elapsed before the key's behaviour started for the first time.
<i>Calculation method</i>	The value of the test clock at the first key press.

<i>Analysis in zones</i>	The value of the test clock at the first key press that occurred within the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

6.4 Latency to first release

<i>Description</i>	Reports the amount of time that elapsed in the test before the key was released for the first time, i.e. the amount of time that elapsed before the key's behaviour ended for the first time.
<i>Calculation method</i>	The value of the test clock at the first key release.
<i>Analysis in zones</i>	The value of the test clock at the first key release that occurred within the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

6.5 Distance travelled before 1st press

<i>Description</i>	Reports the distance the animal had travelled in the apparatus up to the moment the key was first pressed.
<i>Calculation method</i>	The accumulated total distance travelled is noted at the time of the first key press.
<i>Analysis in zones</i>	The measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	None

6.6 Longest press

<i>Description</i>	Reports the duration of the longest period for which the key was continuously pressed down, i.e. the longest period for which the animal continuously exhibited the key's behaviour.
<i>Calculation method</i>	The duration of each key press is calculated and the largest value is found.
<i>Analysis in zones</i>	The longest period for which the key was continuously pressed in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the longest period for which the key was continuously pressed during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

6.7 Shortest press

<i>Description</i>	Reports the duration of the shortest period for which the key was continuously pressed down, i.e. the shortest period for which the animal continuously exhibited the key's behaviour.
<i>Calculation method</i>	The duration of each key press is calculated and the smallest value is found.
<i>Analysis in zones</i>	The shortest period for which the key was continuously pressed in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest period for which the key was continuously pressed during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

6.8 Average press duration

<i>Description</i>	Reports the average duration for which the key was held down, i.e. the average duration of the episodes of the key's behaviour.
<i>Calculation method</i>	Calculated by dividing the Time pressed by the Number of presses.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

6.9 Press frequency

<i>Description</i>	Reports the frequency with which the key was pressed, i.e. the frequency with which the animal exhibited the key's behaviour.
<i>Calculation method</i>	Calculated by dividing the Number of presses by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of presses in the zone by the <i>Total time in the zone</i> . For an investigation zone, this will be the number of presses while the animal was investigating the zone divided by the total time investigating the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of presses which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

6.10 List of activation durations

<i>Description</i>	Reports a comma-separated list of the duration of each activation of the key. For example, if key was activated three times, for 10s on the first occasion and for 20s on the second and third, then the list would show '10.0, 20.0, 20.0'.
<i>Calculation method</i>	The duration of each activation is calculated and added to the list. If the key is active at the end of the test, then the last duration in the list will be from the time the key was activated to the time of the test end.
<i>Analysis in zones</i>	The measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	All the durations in the list are reported in seconds.
<i>Notes</i>	<p>When included on the Data page, this measure will show all the activation durations in a single cell. If the spreadsheet is saved in CSV format and then opened in (for example) Microsoft Excel, then the durations will be listed in individual cells.</p> <p>The length of the list is limited to 8192 characters, but usually at last 1,000 activations will be listed before the limit is reached.</p>

7 On/off input measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the input activation occurred) and also for each defined zone.

7.1 Number of activations

<i>Description</i>	Reports the number of times the input was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	<i>Activation</i> is defined when you create the input; see Setting up an on/off input.

7.2 Time active

<i>Description</i>	Reports the total amount of time that the input was active.
<i>Calculation method</i>	Sums the amount of time for which the input was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the input was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the input is active. In this case, the time the input is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the input was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the input is already active at the start of the period. In this case, the time the input is active will be registered but the input activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	<i>Activation</i> is defined when you create the input; see Setting up an on/off input.

7.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the input was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the input.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the input that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time segment or time period is the latency to the first activation within the time segment or period. If the

input is not activated within the time segment or period then the result is #N/A unless the analysis option to *Use the test duration as the latency for events which don't occur* is set, in which case the result will be the duration of the time segment or period.

Note that an input which is already active at the start of a time segment or time period will NOT be reported as having a latency of zero; this measure reports the latency to the first activation DURING the time segment or period.

Units

Seconds

Notes

Activation is defined when you create the input; see Setting up an on/off input.

7.4 Latency to first deactivation

Description

Reports the amount of time that elapsed in the test before the input was deactivated for the first time.

Calculation method

The value of the test clock at the first input deactivation.

Analysis in zones

The value of the test clock at the first input deactivation that occurred while the animal was within the zone.

Analysis across time

This measure can be analysed across time. The result for a time segment or time period is the latency to the first deactivation within the time segment or period. If the input is not deactivated within the time segment or period then the result is #N/A unless the analysis option to *Use the test duration as the latency for events which don't occur* is set, in which case the result will be the duration of the time segment or period.

Units

Seconds

Notes

Activation is defined when you create the input; see Setting up an on/off input.

7.5 Longest activation

Description

Reports the duration of the longest period for which the input was continuously active.

Calculation method

The duration of each input activation is calculated and the largest value is found.

Analysis in zones

The longest period for which the input was continuously active while the animal was in the zone.

Analysis across time

This measure can be analysed across time. For any time period, the result is the longest period for which the input was continuously active during the period.

Units

Seconds

Notes

Activation is defined when you create the input; see Setting up an on/off input.

7.6 Shortest activation

Description

Reports the duration of the shortest period for which the input was continuously active.

<i>Calculation method</i>	The duration of each input activation is calculated and the smallest value is found.
<i>Analysis in zones</i>	The shortest period for which the input was continuously active while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest period for which the input was continuously active during the period.
<i>Units</i>	Seconds
<i>Notes</i>	<i>Activation</i> is defined when you create the input; see Setting up an on/off input.

7.7 Average activation duration

<i>Description</i>	Reports the average duration for which the input was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	<i>Activation</i> is defined when you create the input; see Setting up an on/off input.

7.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the input was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	<i>Activation</i> is defined when you create the input; see Setting up an on/off input.

8 Rotary encoder measures

Each rotary encoder measure is available for the apparatus as whole, and also for each defined zone.

8.1 Number of rotations

<i>Description</i>	Reports the number of complete rotations of the encoder.
<i>Calculation method</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' same direction pulses is received.
<i>Analysis in zones</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' same direction pulses is received and the animal is in the zone for the entire sequence.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of rotations that ENDED during the time period.
<i>Units</i>	None
<i>Notes</i>	None

8.2 Time turning

<i>Description</i>	Reports the time for which the encoder was actually turning.
<i>Calculation method</i>	Reports the time for which the instantaneous rotational velocity (IRV) was not zero. The IRV is calculated as follows: Starts with a single pulse of a certain direction. Then counts the number of consecutive same direction pulses until at least 200 milliseconds have elapsed. Uses the elapsed time and the number of pulses detected during it to calculate the rotational velocity. This value is averaged (using a moving average) over ten values, and this is the instantaneous rotational velocity.
<i>Analysis in zones</i>	The time the encoder was turning while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the time that the encoder was turning during the time period.
<i>Units</i>	Seconds
<i>Notes</i>	None

8.3 Distance (only applies to running wheels)

<i>Description</i>	Reports the distance the animal 'travelled'.
<i>Calculation method</i>	The number of rotations multiplied by the circumference of the wheel.
<i>Analysis in zones</i>	The distance the animal 'travelled' in the wheel while it was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of rotations that occurred during the time period multiplied by the circumference of the wheel.
<i>Units</i>	Metres
<i>Notes</i>	This result is only available if the device which includes the rotary encoder reports a circumference - this feature is only supported by some running wheels.

8.4 Number of clockwise rotations

<i>Description</i>	Reports the number of complete clockwise rotations of the encoder.
<i>Calculation method</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' clockwise pulses is received.
<i>Analysis in zones</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' clockwise pulses is received and the animal is in the zone for the entire sequence.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of clockwise rotations that ENDED during the time period.
<i>Units</i>	None
<i>Notes</i>	None

8.5 Number of anti-clockwise rotations

<i>Description</i>	Reports the number of complete anti-clockwise rotations of the encoder.
<i>Calculation method</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' anti-clockwise pulses is received.
<i>Analysis in zones</i>	Counts a rotation when an unbroken sequence of 'number of pulses per rotation' anti-clockwise pulses is received and the animal is in the zone for the entire sequence.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of anti-clockwise rotations that ENDED during the time period.
<i>Units</i>	None
<i>Notes</i>	None

8.6 Number of reversals

<i>Description</i>	Reports the number of times the direction of the encoder changed.
<i>Calculation method</i>	Counts the number of times a clockwise pulse was followed by an anti-clockwise pulse and vice versa.
<i>Analysis in zones</i>	Counts the number of times a clockwise pulse was followed by an anti-clockwise pulse and vice versa when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of reversals that occurred during the time period.
<i>Units</i>	None
<i>Notes</i>	None

8.7 Number of half rotations

<i>Description</i>	Reports the number of half rotations of the encoder.
--------------------	--

<i>Calculation method</i>	Counts a half rotation when an unbroken sequence of 'number of pulses per rotation' divided by two, same direction pulses is received.
<i>Analysis in zones</i>	Counts a half rotation when an unbroken sequence of 'number of pulses per rotation' divided by two, same direction pulses is received and the animal is in the zone for the entire sequence.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of half rotations that ENDED during the time period.
<i>Units</i>	None
<i>Notes</i>	If the number of pulses per rotation is not exactly divisible by 2, then the result is rounded down. This can yield some inaccuracy. Consider when the number of pulses per rotation is 5, then the number of pulses per half rotation will be considered to be $5/2 = 2.5$, which when rounded down = 2. So after 2 same direction pulses, a half rotation will be counted; thus after 2 complete rotations 5 half rotations will be counted, when the correct value is 4. For this reason, it is not recommended to use encoders with a number of pulses per rotation (PPR) that is not divisible by 2. Fortunately, many encoders use binary powers for their PPR, for example 16, 32, 64, etc.

8.8 Number of quarter rotations

<i>Description</i>	Reports the number of quarter rotations of the encoder.
<i>Calculation method</i>	Counts a quarter rotation when an unbroken sequence of 'number of pulses per rotation' divided by four, same direction pulses is received.
<i>Analysis in zones</i>	Counts a quarter rotation when an unbroken sequence of 'number of pulses per rotation' divided by four, same direction pulses is received and the animal is in the zone for the entire sequence.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of quarter rotations that ENDED during the time period.
<i>Units</i>	None
<i>Notes</i>	If the number of pulses per rotation is not exactly divisible by 4 then the result is rounded down. This can yield some inaccuracy. Consider when the number of pulses per rotation is 5, then the number of pulses per quarter rotation will be considered to be $5/4 = 1.25$, which when rounded down = 1. So after 1 pulse, a quarter rotation will be counted; thus after 4 complete rotations 20 quarter rotations will be counted, when the correct value is 16. For this reason, it is not recommended to use encoders with a number of pulses per rotation that is not divisible by 4. Fortunately, many encoders use binary powers for their Pulses Per Rotation, for example 16, 32, 64, etc.

8.9 Degrees of clockwise rotation

<i>Description</i>	Reports the number of degrees of clockwise rotation of the encoder.
<i>Calculation method</i>	For each clockwise pulse of the encoder, adds '360 / Number of Pulses Per Rotation' to the result.

<i>Analysis in zones</i>	The number of degrees of clockwise rotation of the encoder while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is updated for each clockwise pulse that occurred during the time period.
<i>Units</i>	Degrees
<i>Notes</i>	This result is useful because it provides fine resolution (compared to number of rotations), while being normalised for all encoders (unlike number of pulses).

8.10 Degrees of anti-clockwise rotation

<i>Description</i>	Reports the number of degrees of anti-clockwise rotation of the encoder.
<i>Calculation method</i>	For each anti-clockwise pulse of the encoder, adds '360 / Number of Pulses Per Rotation' to the result.
<i>Analysis in zones</i>	The number of degrees of anti-clockwise rotation of the encoder while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is updated for each anti-clockwise pulse that occurred during the time period.
<i>Units</i>	Degrees
<i>Notes</i>	This result is useful because it provides fine resolution (compared to number of rotations), while being normalised for all encoders (unlike number of pulses).

8.11 Maximum RPM

<i>Description</i>	Reports the maximum rotational velocity of the encoder in units of revolutions per minute.
<i>Calculation method</i>	Starts with a single pulse of a certain direction, then counts the number of consecutive same-direction pulses until at least 200 milliseconds have elapsed. Uses the elapsed time and the number of pulses detected during it to calculate the rotational velocity. This value is averaged (using a moving average) over ten values, and this is the instantaneous rotational velocity. The highest instantaneous rotational velocity is the maximum RPM.
<i>Analysis in zones</i>	The highest instantaneous rotational velocity while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the maximum value determined during the time period.
<i>Units</i>	Revolutions per minute
<i>Notes</i>	A rotary encoder has an inherent maximum rotational velocity, above which the encoder will not be read accurately. This value is reported on the I/O page when an encoder is selected. In fact, strictly speaking, this is a limitation of the interface that reads the encoder and not of the encoder itself. For example, when reading a 32 PPR encoder with the ANY-maze interface, the maximum RPM is 416, which is approximately 7 revolutions per second (enough for most likely situations in behavioural research). Note that using an encoder with half as many pulses per rotation will double this value.

8.12 Minimum RPM

<i>Description</i>	Reports the minimum rotational velocity of the encoder in units of revolutions per minute.
<i>Calculation method</i>	Starts with a single pulse of a certain direction. Then counts the number of consecutive same direction pulses until at least 200 milliseconds have elapsed. Uses the elapsed time and the number of pulses detected during it to calculate the rotational velocity. This value is averaged (using a moving average) over ten values, and this is the instantaneous rotational velocity. The lowest instantaneous rotational velocity is the minimum RPM.
<i>Analysis in zones</i>	The lowest instantaneous rotational velocity while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the minimum value determined during the time period.
<i>Units</i>	Revolutions per minute
<i>Notes</i>	If the encoder stops turning, this value will be zero.

8.13 Average RPM

<i>Description</i>	Reports the average rotational velocity of the encoder in units of revolutions per minute.
<i>Calculation method</i>	Starts with a single pulse of a certain direction. Then counts the number of consecutive same direction pulses until at least 200 milliseconds have elapsed. Uses the elapsed time and the number of pulses detected during it to calculate the rotational velocity. This value is averaged (using a moving average) over ten values, and this is the instantaneous rotational velocity. The instantaneous rotational velocity is averaged throughout the test to yield this measure (see notes).
<i>Analysis in zones</i>	The average of the instantaneous rotational velocity values reported while the animal was in the zone (see notes).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the average of the instantaneous rotational velocity values reported during the time period (see notes).
<i>Units</i>	Revolutions per minute
<i>Notes</i>	The averaging of the instantaneous rotational velocity (IRV) values uses a time-based averaging technique, such that the average is the sum of each IRV value multiplied by the time between two consecutive values divided by the sum of the time between all values. This is required because the IRV values are not reported at a fixed frequency.

8.14 Average RPM while turning

<i>Description</i>	Reports the average rotational velocity of the encoder for the time when it was turning, in units of revolutions per minute.
--------------------	--

<i>Calculation method</i>	Starts with a single pulse of a certain direction. Then counts the number of consecutive same direction pulses until at least 200 milliseconds have elapsed. Uses the elapsed time and the number of pulses detected during it to calculate the rotational velocity. This value is averaged (using a moving average) over ten values, and this is the instantaneous rotational velocity (IRV). All non-zero IRV values are averaged to yield this measure (see notes).
<i>Analysis in zones</i>	The average of the non-zero instantaneous rotational velocity values reported while the animal was in the zone (see notes).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the average of the non-zero instantaneous rotational velocity values reported during the time period (see notes).
<i>Units</i>	Revolutions per minute
<i>Notes</i>	The averaging of the instantaneous rotational velocity (IRV) values uses a time-based averaging technique, such that the average is the sum of each non-zero IRV value multiplied by the time between two consecutive non-zero values divided by the sum of the time between all such values. This is required because the IRV values are not reported at a fixed frequency.

9 Signal measures

9.1 Average

<i>Description</i>	The mean of the signal's value over the duration of the test.
<i>Calculation method</i>	All the signal's samples are summed across the duration of the test and the sum is divided by the count of samples.
<i>Analysis in zones</i>	This measure can be analysed in zones. It reports the mean of the signal's value while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It reports the mean of the signal's value during the time period.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.2 Maximum

<i>Description</i>	The signal's maximum value over the duration of the test.
<i>Calculation method</i>	All the signal's samples are scanned and the largest value is found.
<i>Analysis in zones</i>	This measure can be analysed in zones. It reports the signal's maximum value reports while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It reports the signal's maximum value during the time period.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.3 Minimum

<i>Description</i>	The signal's minimum value over the duration of the test.
<i>Calculation method</i>	All the signal's samples are scanned and the smallest value is found.
<i>Analysis in zones</i>	This measure can be analysed in zones. It reports the signal's minimum value reports while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It reports the signal's minimum value during the time period.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.4 Time of maximum

<i>Description</i>	The time at which the signal's maximum value occurred
<i>Calculation method</i>	The test clock at the time the maximum value occurred

<i>Analysis in zones</i>	This measure can be analysed in zones. It reports the time at which the signal's maximum value occurred, while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It reports the time <i>in the time period</i> at which the signal's maximum value during the time period occurred.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.5 Time of minimum

<i>Description</i>	The time at which the signal's minimum value occurred
<i>Calculation method</i>	The test clock at the time the minimum value occurred
<i>Analysis in zones</i>	This measure can be analysed in zones. It reports the time at which the signal's minimum value occurred, while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It reports the time <i>in the time period</i> at which the signal's minimum value during the time period occurred.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.6 Baseline

<i>Description</i>	The signal's baseline value
<i>Calculation method</i>	The average of the signal's values during the baseline period. All the signal's samples are summed across the baseline period and the sum is divided by the count of samples.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.7 Baseline deviation

<i>Description</i>	The signal's baseline deviation value, as defined in the protocol. This is the amount the signal must deviate from the baseline value for it to count as a <i>deviation</i> . For example, if the protocol specifies that the signal's baseline deviation is 10% and the baseline value is 20, then the baseline deviation will be 2.
<i>Calculation method</i>	Depends on how the deviation is defined in the protocol.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.8 Baseline standard deviation

<i>Description</i>	The standard deviation of the signal's value during the baseline period.
<i>Calculation method</i>	The standard deviation of all of the signal's samples during the baseline period.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.9 Time of the end of the baseline period

<i>Description</i>	The time at which the baseline period ended
<i>Calculation method</i>	The test clock at the time the baseline period ended
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	If the baseline period is defined in the protocol as a set period, then this will be that period, but if the baseline period ends at some time marker (which could be different in different tests), then this will be the time of the time marker.

9.10 Time of first positive deviation from baseline

<i>Description</i>	The time at which the signal's value first deviated above the baseline value by more than the Baseline deviation.
<i>Calculation method</i>	The test clock at the time at which the deviation occurred.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure can be analysed across time. It is the time <i>within the period</i> at which the signal's value first (within the period) deviated above the baseline value by more than the baseline deviation.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.11 Time of first return to baseline from positive deviation

<i>Description</i>	The time at which the signal's value first returned from being above the baseline value to be within the baseline deviation of it.
<i>Calculation method</i>	The test clock at the time at which the return to the baseline occurred. It is implicit that a positive deviation from the baseline must already have occurred.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.

<i>Analysis across time</i>	This measure can be analysed across time. It is the time <i>within the period</i> at which the signal's value first (within the period) returned from being above the baseline value to be within the baseline deviation of it.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.12 Time of first negative deviation from baseline

<i>Description</i>	The time at which the signal's value first deviated below the baseline value by more than the baseline deviation.
<i>Calculation method</i>	The test clock at the time at which the deviation occurred.
<i>Analysis in zones</i>	This value measure be analysed in zones.
<i>Analysis across time</i>	This measure can be analysed across time. It is the time <i>within the period</i> at which the signal's value first (within the period) deviated below the baseline value by more than the baseline deviation.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.13 Time of first return to baseline from negative deviation

<i>Description</i>	The time at which the signal's value first returned from being below the baseline value to be within the baseline deviation of it.
<i>Calculation method</i>	The test clock at the time at which the return to the baseline occurred. It is implicit that a negative deviation from the baseline must already have occurred.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure can be analysed across time. It is the time <i>within the period</i> at which the signal's value first (within the period) returned from being below the baseline value to be within the baseline deviation of it.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.14 Integral above baseline

<i>Description</i>	The sum of the area under the signal's line and above the baseline
<i>Calculation method</i>	Every sample after the baseline period is compared to the baseline; if it is above it, then the difference between the signal value and the baseline (divided by the sampling rate) is added to the integral.
<i>Analysis in zones</i>	This measure can be analysed in zones. It is the sum of the area under the signal's line and above the baseline for all the samples captured while the animal was in the zone.

<i>Analysis across time</i>	This measure can be analysed across time. It is the sum of the area under the signal's line and above the baseline for all the samples captured during the time period.
<i>Units</i>	Same as the signal's units · seconds
<i>Notes</i>	None

9.15 Integral below baseline

<i>Description</i>	The sum of the area above the signal's line and below the baseline
<i>Calculation method</i>	Every sample after the baseline period is compared to the baseline; if it is below it, then the difference between the baseline and the signal value (divided by the sampling rate) is added to the integral.
<i>Analysis in zones</i>	This measure can be analysed in zones. It is the sum of the area above the signal's line and below the baseline for all the samples captured while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. It is the sum of the area above the signal's line and below the baseline for all the samples captured during the time period.
<i>Units</i>	Same as the signal's units · seconds
<i>Notes</i>	None

9.16 Average maximum for each visit to zone

<i>Description</i>	The average of the signal's maximum value on each of the animal's visits to the zone. This measure can only be reported for zones - it cannot be reported for the apparatus as a whole.
<i>Calculation method</i>	The maximum signal value on each of the animal's visits to the zone is found. These are summed and divided by the number of zone visits.
<i>Analysis in zones</i>	This measure can only be reported in zones.
<i>Analysis across time</i>	This measure can be analysed across time. For each time period, it is the average of the signal's maximum value on each of the animal's visits to the zone during the period.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.17 Average time to maximum for each visit to zone

<i>Description</i>	The average time after entering a zone at which the signal reaches the maximum value during the visit to the zone. This measure can only be reported for zones - it cannot be reported for the apparatus as a whole.
<i>Calculation method</i>	For each visit to the zone, the maximum signal value is found. The value of the test clock when the animal entered the zone is subtracted from the value of the test

clock at the time of this maximum, and the resulting time is added to a sum of times. This sum is then divided by the total number of zone visits.

<i>Analysis in zones</i>	This measure can only be reported in zones.
<i>Analysis across time</i>	This measure can be analysed across time. For each visit to the zone during the time period, the maximum signal value is found. The value of the test clock when the animal entered the zone (or when the time period started, whichever is later) is subtracted from the value of the test clock at the time of this maximum, and the resulting time is added to a sum of times. This sum is then divided by the total number of zone visits during the time period.
<i>Units</i>	Seconds
<i>Notes</i>	None

9.18 Average minimum for each visit to zone

<i>Description</i>	The average of the signal's minimum value on each of the animal's visits to the zone. This measure can only be reported for zones - it cannot be reported for the apparatus as a whole.
<i>Calculation method</i>	The minimum signal value on each of the animal's visits to the zone is found. These are summed and divided by the number of zone visits.
<i>Analysis in zones</i>	This measure can only be reported in zones.
<i>Analysis across time</i>	This measure can be analysed across time. For each time period, it is the average of the signal's minimum value on each of the animal's visits to the zone during the period.
<i>Units</i>	Same as the signal's units
<i>Notes</i>	None

9.19 Average time to minimum for each visit to zone

<i>Description</i>	The average time after entering a zone at which the signal reaches the minimum value during the visit to the zone. This measure can only be reported for zones - it cannot be reported for the apparatus as a whole.
<i>Calculation method</i>	For each visit to the zone, the minimum signal value is found. The value of the test clock when the animal entered the zone is subtracted from the value of the test clock at the time of this minimum, and the resulting time is added to a sum of times. This sum is then divided by the total number of zone visits.
<i>Analysis in zones</i>	This measure can only be reported in zones.
<i>Analysis across time</i>	This measure can be analysed across time. For each visit to the zone during a time period, the minimum signal value is found. The value of the test clock when the animal entered the zone (or when the time period started, whichever is later) is subtracted from the value of the test clock at the time of this minimum, and the resulting time is added to a sum of times. This sum is then divided by the total number of zone visits during the time period.
<i>Units</i>	Seconds

Notes None

9.20 Average at zone entry

Description The average of the signal's values at the moment the animal entered the zone. This measure can **only** be reported for zones - it cannot be reported for the apparatus as a whole.

Calculation method Each time the animal enters the zone, the signal's value is added to a sum. This sum is then divided by the number of zone entries.

Analysis in zones This measure can **only** be reported in zones.

Analysis across time This measure can be analysed across time. For each zone entry during a time period, the signal's value at the moment the animal entered is added to a sum. This sum is then divided by the number of zone entries during the period.

Units Same as the signal's units

Notes None

9.21 Average at zone exit

Description The average of the signal's values at the moment the animal exited the zone. This measure can **only** be reported for zones - it cannot be reported for the apparatus as a whole.

Calculation method Each time the animal leaves the zone, the signal's value is added to a sum. This sum is then divided by the number of zone exists.

Analysis in zones This measure can **only** be reported in zones.

Analysis across time This measure can be analysed across time. For each zone exit during a time period, the signal's value at the moment the animal exited is added to a sum. This sum is then divided by the number of zone exists during the period.

Units Same as the signal's units

Notes None

10 Sensor measures

10.1 Initial value

<i>Description</i>	The sensor's value at the start of the test.
<i>Calculation method</i>	Depends on how the sensor is configured; either the last value reported by the sensor prior to the start of the test (this will typically be read less than 1 second before the start of the test), or the value read from the sensor when the sensor's key was pressed before the test started.
<i>Analysis in zones</i>	This value is not available within zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Depends on the sensor (lux, degrees, grams or %)
<i>Notes</i>	This is the only measure that is reported when a sensor is set to be read at the test start. When a sensor is set to report the change in its value during a test, this measure is NOT reported (as it would always be 0).

10.2 Average value

<i>Description</i>	Reports the average value reported by the sensor.
<i>Calculation method</i>	Simple average of the values reported by the sensor during the test.
<i>Analysis in zones</i>	Average of the values reported while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the average of the values reported during the time period.
<i>Units</i>	Depends on the sensor (lux, degrees, grams or %)
<i>Notes</i>	This measure is only reported when a sensor is set to be read continuously throughout a test.

10.3 Maximum value

<i>Description</i>	Reports the maximum value reported by the sensor.
<i>Calculation method</i>	The largest value reported by the sensor during the test.
<i>Analysis in zones</i>	The largest value reported while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the largest of the values reported during the time period.
<i>Units</i>	Depends on the sensor (lux, degrees, grams or %)
<i>Notes</i>	This measure is only reported when a sensor is set to be read continuously throughout a test.

10.4 Minimum value

<i>Description</i>	Reports the minimum value reported by the sensor.
--------------------	---

<i>Calculation method</i>	The smallest value reported by the sensor during the test.
<i>Analysis in zones</i>	The smallest value reported while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the smallest of the values reported during the time period.
<i>Units</i>	Depends on the sensor (lux, degrees, grams or %)
<i>Notes</i>	This measure is only reported when a sensor is set to be read continuously throughout a test.

10.5 Change

<i>Description</i>	Reports the change in the sensor's value relative to the value at the start of the test.
<i>Calculation method</i>	The first value recorded in the test is subtracted from the last value recorded in the test.
<i>Analysis in zones</i>	This value is not available within zones.
<i>Analysis across time</i>	This measure can be analysed across time (see notes). The result is the value at the start of the time period subtracted from the last value recorded in the time period.
<i>Units</i>	Depends on the sensor (lux, degrees, grams or %)
<i>Notes</i>	This measure is only reported when a sensor is set to be read continuously throughout a test or is set to be read after the test when the sensor's key is pressed. If the sensor is not set to be read continuously, then analysis across time is not possible.

11 Movement detector measures

11.1 Count

<i>Description</i>	Reports the number of times beams in the movement detector's photobeam array were broken.
<i>Calculation method</i>	Counts breaks reported by the photobeam array. Note that ANY-maze does NOT count repeated breaks of the same beam - this is to avoid counting small oscillatory movements of the animal's body as <i>movement</i> (i.e. as in changes of location). For example, consider an animal that moves across the area covered by a photobeam array: we would see beam 1 broken, then beam 2, then beam 3, etc. Now consider an animal that is sitting in one place and grooming; its leg might keep breaking the same beam again and again; this would not be counted as movement.
<i>Analysis in zones</i>	Not available
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of counts that occurred during the time period.
<i>Units</i>	None
<i>Notes</i>	None

11.2 Time moving

<i>Description</i>	Reports the amount of time that the animal was moving.
<i>Calculation method</i>	Sums the durations of all the individual bouts of movement.
<i>Analysis in zones</i>	Not available
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the duration of the movement within the time period.
<i>Units</i>	Seconds
<i>Notes</i>	Movements are detected by a beam being broken. This is necessarily an instantaneous event, but it is deemed to mean that the animal is moving and will continue to do so for a certain period (the movement detector's time-out). If another beam is broken within this period, then the animal is deemed to still be moving and so the movement bout extends to include this break. This continues until no break occurs within the timeout period, when the bout is then deemed to end. This mechanism is described in detail in this topic.

11.3 Latency to first beam break

<i>Description</i>	Reports the latency to the first beam being broken in the test - i.e. the first bout of movement.
<i>Calculation method</i>	The time up to the first beam break <i>during</i> the test is measured.
<i>Analysis in zones</i>	Not available
<i>Analysis across time</i>	This measure cannot be analysed across time.

Units
Notes

Seconds

Beams that are already broken at the test start have no effect on this measure.

12 On/off output measures

12.1 Number of activations

<i>Description</i>	Reports the number of times the on/off output was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

12.2 Time active

<i>Description</i>	Reports the total amount of time that the on/off output was active.
<i>Calculation method</i>	Sums the amount of time for which the on/off output was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the on/off output was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the output is active. In this case, the time the output is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the output was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the output is already active at the start of the period. In this case, the time the output is active will be registered but the activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

12.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the on/off output was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the on/off output.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the on/off output that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

12.4 Latency to first deactivation

<i>Description</i>	Reports the amount of time that elapsed in the test before the on/off output was deactivated for the first time.
<i>Calculation method</i>	The value of the test clock at the first on/off output deactivation.
<i>Analysis in zones</i>	The value of the test clock at the first on/off output deactivation that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

12.5 Longest activation

<i>Description</i>	Reports the duration of the longest period for which the on/off output was continuously active.
<i>Calculation method</i>	The duration of each output activation is calculated and the largest value is found.
<i>Analysis in zones</i>	The longest period for which the on/off output was continuously active while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the longest period for which the on/off output was continuously active during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

12.6 Shortest activation

<i>Description</i>	Reports the duration of the shortest period for which the on/off output was continuously active.
<i>Calculation method</i>	The duration of each output activation is calculated and the smallest value is found.
<i>Analysis in zones</i>	The shortest period for which the on/off output was continuously active while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest period for which the on/off output was continuously active during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

12.7 Average activation duration

<i>Description</i>	Reports the average duration for which the on/off output was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.

<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

12.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the on/off output was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

13 Speaker measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the speaker activation occurred), and also for each defined zone.

13.1 Number of activations

<i>Description</i>	Reports the number of times the speaker was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

13.2 Time active

<i>Description</i>	Reports the total amount of time that the speaker was active.
<i>Calculation method</i>	Sums the amount of time for which the speaker was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the speaker was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the speaker is active. In this case, the time the speaker is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the speaker was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the speaker is already active at the start of the period. In this case, the time the speaker is active will be registered but the activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

13.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the speaker was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the speaker.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the speaker that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds

Notes None

13.4 Latency to first deactivation

Description Reports the amount of time that elapsed in the test before the speaker was deactivated for the first time.

Calculation method The value of the test clock at the first deactivation of the speaker.

Analysis in zones The value of the test clock at the first deactivation of the speaker that occurred while the animal was within the zone.

Analysis across time This measure cannot be analysed across time.

Units Seconds

Notes None

13.5 Longest activation

Description Reports the duration of the longest period for which the speaker was continuously active.

Calculation method The duration of each speaker activation is calculated and the largest value is found.

Analysis in zones The longest period for which the speaker was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the longest period for which the speaker was continuously active during the period.

Units Seconds

Notes None

13.6 Shortest activation

Description Reports the duration of the shortest period for which the speaker was continuously active.

Calculation method The duration of each speaker activation is calculated and the smallest value is found.

Analysis in zones The shortest period for which the speaker was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the shortest period for which the speaker was continuously active during the period.

Units Seconds

Notes None

13.7 Average activation duration

<i>Description</i>	Reports the average duration for which the speaker was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

13.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the speaker was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

14 Shocker measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the shocker activation occurred), and also for each defined zone.

14.1 Number of activations

<i>Description</i>	Reports the number of times the shocker was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

14.2 Time active

<i>Description</i>	Reports the total amount of time that the shocker was active.
<i>Calculation method</i>	Sums the amount of time for which the shocker was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the shocker was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the shocker is active. In this case, the time the shocker is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the shocker was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the shocker is already active at the start of the period. In this case, the time the shocker is active will be registered but the activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

14.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the shocker was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the shocker.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the shocker that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds

Notes None

14.4 Latency to first deactivation

Description Reports the amount of time that elapsed in the test before the shocker was deactivated for the first time.

Calculation method The value of the test clock at the first deactivation of the shocker.

Analysis in zones The value of the test clock at the first deactivation of the shocker that occurred while the animal was within the zone.

Analysis across time This measure cannot be analysed across time.

Units Seconds

Notes None

14.5 Longest activation

Description Reports the duration of the longest period for which the shocker was continuously active.

Calculation method The duration of each shocker activation is calculated and the largest value is found.

Analysis in zones The longest period for which the shocker was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the longest period for which the shocker was continuously active during the period.

Units Seconds

Notes None

14.6 Shortest activation

Description Reports the duration of the shortest period for which the shocker was continuously active.

Calculation method The duration of each shocker activation is calculated and the smallest value is found.

Analysis in zones The shortest period for which the shocker was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the shortest period for which the shocker was continuously active during the period.

Units Seconds

Notes None

14.7 Average activation duration

<i>Description</i>	Reports the average duration for which the shocker was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

14.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the shocker was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

15 Pellet dispenser measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the pellet dispenser activation occurred), and also for each defined zone.

15.1 Number of pellets dispensed

<i>Description</i>	Reports the number of pellets that were dispensed. For Ugo Basile pellet dispensers, feedback from the device ensures that this is the actual count of pellets dispensed; for all other pellet dispensers, this is the count of the number of times that ANY-maze asked the device to dispense a pellet (since there is no feedback, ANY-maze must assume that this was successful!).
<i>Calculation method</i>	Counts the number of pellets dispensed.
<i>Analysis in zones</i>	Counts the number of pellets dispensed when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of pellets dispensed during the time period.
<i>Units</i>	None
<i>Notes</i>	None

15.2 Latency to first pellet dispensed

<i>Description</i>	Reports the amount of time that elapsed in the test before the first pellet was dispensed.
<i>Calculation method</i>	The value of the test clock when the first pellet was dispensed.
<i>Analysis in zones</i>	The value of the test clock at the time the first pellet was dispensed while the animal was within the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The value is the time since the start of the time segment or time period that the first pellet was dispensed in the time segment or time period. If no pellet was dispensed in the time segment or time period then the result is #N/A, unless the analysis option to <i>Use the test duration as the latency for events which don't occur</i> is set, in which case the result will be the duration of the time segment or period.
<i>Units</i>	Seconds
<i>Notes</i>	None

16 Laser controller measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the laser controller activation occurred), and also for each defined zone.

16.1 Number of activations

<i>Description</i>	Reports the number of times the laser controller was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

16.2 Time active

<i>Description</i>	Reports the total amount of time that the laser controller was active.
<i>Calculation method</i>	Sums the amount of time for which the laser controller was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the laser controller was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the controller is active. In this case, the time the controller is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the laser controller was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the controller is already active at the start of the period. In this case, the time the controller is active will be registered but the activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

16.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the laser controller was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the laser controller.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the laser controller that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds

Notes None

16.4 Latency to first deactivation

Description Reports the amount of time that elapsed in the test before the laser controller was deactivated for the first time.

Calculation method The value of the test clock at the first laser controller deactivation.

Analysis in zones The value of the test clock at the first laser controller deactivation that occurred while the animal was within the zone.

Analysis across time This measure cannot be analysed across time.

Units Seconds

Notes None

16.5 Longest activation

Description Reports the duration of the longest period for which the laser controller was continuously active.

Calculation method The duration of each laser controller activation is calculated, and the largest value is found.

Analysis in zones The longest period for which the laser controller was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the longest period for which the laser controller was continuously active during the period.

Units Seconds

Notes None

16.6 Shortest activation

Description Reports the duration of the shortest period for which the laser controller was continuously active.

Calculation method The duration of each laser controller activation is calculated, and the smallest value is found.

Analysis in zones The shortest period for which the laser controller was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the shortest period for which the laser controller was continuously active during the period.

Units Seconds

Notes None

16.7 Average activation duration

<i>Description</i>	Reports the average duration for which the laser controller was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

16.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the laser controller was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period, divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

17 Syringe pump measures

17.1 Volume infused

<i>Description</i>	The total volume infused by the syringe pump.
<i>Calculation method</i>	The sum of the amounts infused each time the pump was run in the infuse direction during the test. The amount infused is calculated by multiplying the infusion rate by the amount of time the pump was running.
<i>Analysis in zones</i>	The volume infused while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the total volume infused during the time period.
<i>Units</i>	Microlitres
<i>Notes</i>	None.

17.2 Volume withdrawn

<i>Description</i>	The total volume withdrawn by the syringe pump.
<i>Calculation method</i>	The sum of the amounts withdrawn each time the pump was run in the withdraw direction during the test. The amount withdrawn is calculated by multiplying the withdrawal rate by the amount of time the pump was running.
<i>Analysis in zones</i>	The volume withdrawn while the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the total volume withdrawn during the time period.
<i>Units</i>	Microlitres
<i>Notes</i>	Not all syringe pumps can be run in the withdraw direction. For those that can't, the result of this measure will always be zero.

18 Light controller measures

Each measure is available for the apparatus as whole (i.e. irrespective of where the animal was when the light controller activation occurred), and also for each defined zone.

18.1 Number of activations

<i>Description</i>	Reports the number of times the light controller was activated.
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of activations when the animal was in the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

18.2 Time active

<i>Description</i>	Reports the total amount of time that the light controller was active.
<i>Calculation method</i>	Sums the amount of time for which the light controller was active.
<i>Analysis in zones</i>	<p>Sums the amount of time for which the light controller was active while the animal was in the zone.</p> <p>For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the light controller is active. In this case, the time the light controller is active will be registered but the activation itself won't be.</p>
<i>Analysis across time</i>	<p>This measure can be analysed across time. For any time period, the result is the amount of time that the light controller was active during the period.</p> <p>For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the light controller is already active at the start of the period. In this case, the time the light controller is active will be registered but the activation itself won't be.</p>
<i>Units</i>	Seconds
<i>Notes</i>	None

18.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the light controller was activated for the first time.
<i>Calculation method</i>	The value of the test clock at the first activation of the light controller.
<i>Analysis in zones</i>	The value of the test clock at the first activation of the light controller that occurred while the animal was within the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds

Notes None

18.4 Latency to first deactivation

Description Reports the amount of time that elapsed in the test before the light controller was deactivated for the first time.

Calculation method The value of the test clock at the first deactivation of the light controller.

Analysis in zones The value of the test clock at the first deactivation of the light controller that occurred while the animal was within the zone.

Analysis across time This measure cannot be analysed across time.

Units Seconds

Notes None

18.5 Longest activation

Description Reports the duration of the longest period for which the light controller was continuously active.

Calculation method The duration of each light controller activation is calculated and the largest value is found.

Analysis in zones The longest period for which the light controller was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the longest period for which the light controller was continuously active during the period.

Units Seconds

Notes None

18.6 Shortest activation

Description Reports the duration of the shortest period for which the light controller was continuously active.

Calculation method The duration of each light controller activation is calculated and the smallest value is found.

Analysis in zones The shortest period for which the light controller was continuously active while the animal was in the zone.

Analysis across time This measure can be analysed across time. For any time period, the result is the shortest period for which the light controller was continuously active during the period.

Units Seconds

Notes None

18.7 Average activation duration

<i>Description</i>	Reports the average duration for which the light controller was active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

18.8 Frequency of activations

<i>Description</i>	Reports the frequency with which the light controller was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> .
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None

19 OPAD measures

19.1 Temperature when contact broken

<i>Description</i>	The average temperature at which the animal broke contact with the thermal elements.
<i>Calculation method</i>	Each time the animal breaks contact, the temperature of the thermal elements is noted. The result is the average of all the noted values.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the average temperature at which the animal broke contact. for all the breaks which occurred during the time period.
<i>Units</i>	Seconds
<i>Notes</i>	If OPAD is being used with independent control of the left and right thermal elements, then this measure will be reported independently for the left and right sides.

19.2 Number of non-lick contacts

<i>Description</i>	Reports the number of licks that occurred when the animal was not making contact with the thermal elements.
<i>Calculation method</i>	Each time the animal licks, if it is not in contact with the thermal elements, a counter is increased.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	None

19.3 Temperature of interest: time in contact

<i>Description</i>	Reports the amount of time the animal was in contact with the thermal elements when the elements' temperature was within the bounds of the specific temperature of interest.
<i>Calculation method</i>	Each time the animal makes contact and the thermal elements are within the bounds of the temperature of interest (or when the animal is already in contact and the temperature changes to be within the bounds of the temperature of interest), a timer is started. Whenever the animal breaks contact (or the animal is in contact and the temperature changes to be outside the bounds of the temperature of interest), then the timer is stopped. The duration of each timed bout is summed to give the total time in contact.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	If OPAD is being used with independent control of the left and right thermal elements, then this measure will be reported independently for the left and right sides.

19.4 Temperature of interest: number of times contact broken

<i>Description</i>	Reports the number of times that the animal broke contact with the thermal elements when the elements' temperature was within the bounds of the specific temperature of interest.
<i>Calculation method</i>	Each time the animal breaks contact with the thermal elements and the elements are within the bounds of the temperature of interest, a counter is increased.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	If OPAD is being used with independent control of the left and right thermal elements, then this measure will be reported independently for the left and right sides.

19.5 Temperature of interest: number of times contact made

<i>Description</i>	Reports the number of times that the animal made contact with the thermal elements when the elements' temperature was within the bounds of the specific temperature of interest.
<i>Calculation method</i>	Each time the animal makes contact with the thermal elements and the elements are within the bounds of the temperature of interest, a counter is increased.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	If OPAD is being used with independent control of the left and right thermal elements, then this measure will be reported independently for the left and right sides.

19.6 Temperature of interest: number of licks

<i>Description</i>	Reports the number of licks that occurred while the thermal elements' temperature was within the bounds of the specific temperature of interest.
<i>Calculation method</i>	Each time the animal licks and the thermal elements are within the bounds of the temperature of interest, a counter is increased.
<i>Analysis across time</i>	This measure can be analysed across time.
<i>Units</i>	None
<i>Notes</i>	If OPAD is being used with independent control of the left and right thermal elements, then this measure will be reported independently for the left and right sides.

20 Procedure measures

Procedures don't have any measures defined themselves, but you can create your own measures by using *result variables*. These are essentially numeric values that your procedure can set as it runs, and these values can then be included in ANY-maze result analysis.

20.1 Value

<i>Description</i>	Reports the value of the result variable at the end of the test.
<i>Calculation method</i>	The value of the result variable at the end of the test. If the value is never noted, the value will be reported as 0 (since result variables are initialised to 0 at the start of a test).
<i>Analysis in zones</i>	This measure cannot be analysed between zones. If the animal is not in the zone at the end of the test, the value will be undefined.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	No units defined; this is dependent on how the procedure calculates the value of this result variable.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'only once, at the end of the test'.

20.2 Average value

<i>Description</i>	Reports the average value of the result variable.
<i>Calculation method</i>	Simple average of the values noted by the procedure during the test. If the value is never noted, the average will be undefined.
<i>Analysis in zones</i>	Average of the values reported while the animal is in the zone. If the value is never noted while the animal is in the zone, the average will be undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the average of the values reported during the time period. If the value is never noted during the time period, the average will be undefined.
<i>Units</i>	No units defined; this is dependent on how the procedure calculates the value of this result variable.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

20.3 Maximum value

<i>Description</i>	Reports the maximum value of the result variable noted by the procedure during the test.
<i>Calculation method</i>	The largest value of the result variable noted by the procedure during the test. If the value is never noted, the maximum value will be 0 (since result variables are initialised to 0 at the start of a test).

<i>Analysis in zones</i>	The largest value reported while the animal is in the zone. If the value is never noted while the animal is in the zone, the maximum value will be undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the largest of the values reported during the time period. If the value is never noted during the time period, the maximum value will be undefined.
<i>Units</i>	No units defined; this is dependent on how the procedure calculates the value of this result variable.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

20.4 Minimum value

<i>Description</i>	Reports the minimum value of the result variable noted by the procedure during the test.
<i>Calculation method</i>	The smallest value of the result variable noted by the procedure during the test. If the value is never noted, the minimum value will be 0 (since result variables are initialised to 0 at the start of a test).
<i>Analysis in zones</i>	The smallest value reported while the animal is in the zone. If the value is never noted while the animal is in the zone, the minimum value will be undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result is the smallest of the values reported during the time period. If the value is never noted during the time period, the minimum value will be undefined.
<i>Units</i>	No units defined; this is dependent on how the procedure calculates the value of this result variable.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

20.5 Sum of values

<i>Description</i>	Reports the total sum of the result variable's values during the test.
<i>Calculation method</i>	Every time the result variable is noted, the value is added to the ongoing sum. If the value is never noted, the sum of values will be 0 (since result variables are initialised to 0 at the start of a test).
<i>Analysis in zones</i>	The sum of the values reported while the animal is in the zone. If the value is never noted while the animal is in the zone, the sum will be 0.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a time period is the total sum of all values of the result variable that are noted during the time period.

If the value is never noted during the time period, the sum of values will be 0 (since a user-defined variable is always initialised to 0).

<i>Units</i>	No units defined; this is dependent on how the procedure calculates the value of this result variable.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

20.6 Count of values

<i>Description</i>	Reports the total number of times during the test that the result variable was noted.
<i>Calculation method</i>	Each time the result variable is noted, a counter is incremented. If the value is never noted, the count of values will be 0.
<i>Analysis in zones</i>	Counts the number of times the result variable is noted while the animal is in the zone. If the value is never noted while the animal is in the zone, the count will be 0.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the number of times the result variable is noted during the period. If the value is never noted during this time period, the count of values will be 0.
<i>Units</i>	None; this is simply a count of the number of times that a result variable was set by a procedure.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

20.7 List of values

<i>Description</i>	Lists all the values of the result variable stored during the test.
<i>Calculation method</i>	Not applicable
<i>Analysis in zones</i>	Lists all the values of the result variable that were stored while the animal was in the zone.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Can be specified as part of the variable's definition.
<i>Notes</i>	This measure is only available if the result variable is set up to be noted 'every time it is set by the procedure' or 'only when explicitly set'.

21 Event measures

Events and actions have been replaced by procedures, but they are documented here as they can still be used in legacy experiments.

21.1 Number of events

<i>Description</i>	Reports the number of times an event occurred.
<i>Calculation method</i>	Each time the event occurs, a counter is incremented.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the number of times the event occurred during the period.
<i>Units</i>	None
<i>Notes</i>	Event measures only report events which were detected while a test was being performed. Unlike almost all other measures, an event measure will not update to reflect post-test changes in the event definition.

21.2 Latency to first event

<i>Description</i>	Reports the latency to the first occurrence of an event.
<i>Calculation method</i>	When the event occurs, the test time is noted as the latency. If the event doesn't occur during the test, then the result is undefined.
<i>Analysis across time</i>	This measure can be analysed across time. The result for a period is the time within the period that the event first occurred. If the event doesn't occur during the time period, then the result is undefined.
<i>Units</i>	Seconds
<i>Notes</i>	Event measures only report events which were detected while a test was being performed. Unlike almost all other measures, an event measure will not update to reflect post-test changes in the event definition.

22 Virtual switch measures

Most Virtual switch measures are available for the apparatus as whole (i.e. irrespective of where the behaviour occurred), and also for each defined zone. The exceptions are 'Average activation duration' and 'Distance travelled before 1st activation', which are only available for the apparatus as a whole.

22.1 Number of activations

<i>Description</i>	Reports the number of times the virtual switch was activated (turned on).
<i>Calculation method</i>	Counts the number of activations.
<i>Analysis in zones</i>	Counts the number of times the virtual switch was activated while the animal was in the zone (or for an investigation zone, while the animal was investigating the zone). Note that if the virtual switch is active when the animal enters (or starts investigating) the zone, an activation will not be counted.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the number of activations during the time period.
<i>Units</i>	None
<i>Notes</i>	None

22.2 Time active

<i>Description</i>	Reports the total amount of time the virtual switch was active.
<i>Calculation method</i>	Sums the duration of each activation of the virtual switch.
<i>Analysis in zones</i>	Sums the amount of time for which the virtual switch was active when the animal was in the zone (or for an investigation zone, while the animal was investigating the zone). For a particular zone, it's possible for the <i>Time active</i> to be non-zero while the Number of activations is zero. This can occur if the animal enters the zone when the virtual switch is active. In this case, the time the virtual switch is active in the zone will be registered, but the activation itself won't be (as it didn't occur in the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the amount of time that the virtual switch was active during the period. For a particular time period, it's possible for the <i>Time active</i> during the period to be non-zero while the Number of activations for the period is zero. This can occur if the virtual switch is already active at the start of the period. In this case the time the virtual switch is active will be registered, but the activation itself won't be.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.3 Latency to first activation

<i>Description</i>	Reports the amount of time that elapsed in the test before the virtual switch was activated for the first time.
--------------------	---

<i>Calculation method</i>	The value of the test clock at the first activation.
<i>Analysis in zones</i>	The value of the test clock at the first activation that occurred when the animal was in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.4 Latency to first deactivation

<i>Description</i>	Reports the amount of time that elapsed in the test before the virtual switch was deactivated for the first time.
<i>Calculation method</i>	The value of the test clock at the first deactivation.
<i>Analysis in zones</i>	The value of the test clock at the first deactivation that occurred when the animal was in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.5 Distance travelled before 1st activation

<i>Description</i>	Reports the distance the animal had travelled in the apparatus up to the moment the virtual switch was first activated.
<i>Calculation method</i>	The accumulated total distance travelled is noted at the time of the first activation.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Metres
<i>Notes</i>	None

22.6 Distance travelled while active

<i>Description</i>	Reports the distance the animal travelled in the apparatus while the virtual switch was active.
<i>Calculation method</i>	While virtual switch is active, the distance travelled is accumulated.
<i>Analysis in zones</i>	Calculated by summing the distance the animal travelled while in the zone (or, for an investigation zone, while investigating the zone) when the virtual switch was active.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the distance the animal travelled during the time period while the virtual switch was active.

<i>Units</i>	Metres
<i>Notes</i>	None

22.7 Longest activation

<i>Description</i>	Reports the duration of the longest period for which the virtual switch was continuously active.
<i>Calculation method</i>	The duration of each activation is calculated and the largest value is found.
<i>Analysis in zones</i>	The longest period for which the virtual switch was continuously active while the animal was in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the longest period for which the virtual switch was continuously active during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.8 Shortest activation

<i>Description</i>	Reports the duration of the shortest period for which the virtual switch was continuously active.
<i>Calculation method</i>	The duration of each activation is calculated and the smallest value is found.
<i>Analysis in zones</i>	The shortest period for which the virtual switch was continuously active while the animal was in the zone (or for an investigation zone, while the animal was investigating the zone).
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the shortest period for which the virtual switch was continuously active during the period.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.9 Average activation duration

<i>Description</i>	Reports the average duration for which the virtual switch was held active.
<i>Calculation method</i>	Calculated by dividing the Time active by the Number of activations.
<i>Analysis in zones</i>	This measure cannot be analysed in zones.
<i>Analysis across time</i>	This measure cannot be analysed across time.
<i>Units</i>	Seconds
<i>Notes</i>	None

22.10 Frequency of activations

<i>Description</i>	Reports the frequency with which the virtual switch was activated.
<i>Calculation method</i>	Calculated by dividing the Number of activations by the <i>Test duration</i> .
<i>Analysis in zones</i>	Calculated by dividing the Number of activations in the zone by the <i>Total time in the zone</i> . For an investigation zone, this will be the number of activations while the animal was investigating the zone divided by the total time investigating the zone.
<i>Analysis across time</i>	This measure can be analysed across time. For any time period, the result is the Number of activations which occurred during the period divided by the period's duration.
<i>Units</i>	Hertz
<i>Notes</i>	None